

## Azlon® Wide Neck Venting Wash Bottles



- Wide neck
- Square shoulders
- DripLok® vapour venting valve which prevents drips
- Colour coded one piece PP closures that maintain its shape in use
- Closures feature precision jet tip aperture for fine control
- LDPE
- Capacity: 500mL

Supplied in packs of 5

Code	Alt Ref	Closure Colour	Price
<b>BOT9150</b>	WGW234PVENT	Orange	£21.70
<b>BOT9152</b>	WGW235PVENT	Blue	£21.70
<b>BOT9154</b>	WGW236PVENT	White	£21.70
<b>BOT9156</b>	WGW237PVENT	Red	£21.70
<b>BOT9158</b>	WGW238PVENT	Yellow	£21.70
<b>BOT9160</b>	WGW239PVENT	Green	£21.70



## Azlon® Multilingual Wash Bottles



- Wide mouth
- Indelibly printed in four languages (English, French, German and Spanish)
- Choice of three water and five solvent labels, all with Globally Harmonised Symbols, NFPA diamond, risk phrases and safety information
- Available with or without DripLok® vapour venting technology that prevents drips
- Colour coded one piece PP closures that maintain their shape in use
- Closures feature precision jet tip aperture for fine control
- LDPE

Supplied in packs of 5

Code	Alt Ref	Capacity, mL	Label	Description	Cap Colour	Price
<b>BOT9330</b>	WGW531VTML	250	Acetone	Venting	Red	£27.50
<b>BOT9356</b>	WGW531PML	250	Acetone	Non venting	Red	£23.50
<b>BOT9332</b>	WGW535VTML	250	Ethanol	Venting	Orange	£27.40
<b>BOT9364</b>	WGW535PML	250	Ethanol	Non venting	Orange	£24.00
<b>BOT9334</b>	WGW534VTML	250	I.M.S.	Venting	Yellow	£26.70
<b>BOT9362</b>	WGW534PML	250	I.M.S.	Non venting	Yellow	£23.50
<b>BOT9336</b>	WGW533VTML	250	Isopropanol	Venting	Blue	£27.40
<b>BOT9360</b>	WGW533PML	250	Isopropanol	Non venting	Blue	£23.50
<b>BOT9338</b>	WGW532VTML	250	Methanol	Venting	Green	£28.00
<b>BOT9358</b>	WGW532PML	250	Methanol	Non venting	Green	£23.50
<b>BOT9350</b>	WGW600PML	250	Distilled water	Non venting	White	£24.40
<b>BOT9352</b>	WGW602PML	250	Deionised water	Non venting	White	£24.40
<b>BOT9354</b>	WGW530PML	250	Water	Non venting	White	£23.50
<b>BOT9340</b>	WGW537VTML	500	Acetone	Venting	Red	£30.50
<b>BOT9372</b>	WGW537PML	500	Acetone	Non venting	Red	£29.00
<b>BOT9342</b>	WGW541VTML	500	Ethanol	Venting	Orange	£29.20
<b>BOT9380</b>	WGW541PML	500	Ethanol	Non venting	Orange	£26.00
<b>BOT9344</b>	WGW540VTML	500	I.M.S.	Venting	Yellow	£29.70
<b>BOT9378</b>	WGW540PML	500	I.M.S.	Non venting	Yellow	£25.40
<b>BOT9346</b>	WGW539VTML	500	Isopropanol	Venting	Blue	£30.50
<b>BOT9376</b>	WGW539PML	500	Isopropanol	Non venting	Blue	£26.00
<b>BOT9348</b>	WGW538VTML	500	Methanol	Venting	Green	£29.00
<b>BOT9374</b>	WGW538PML	500	Methanol	Non venting	Green	£26.50
<b>BOT9368</b>	WGW601PML	500	Distilled water	Non venting	White	£26.00
<b>BOT9370</b>	WGW603PML	500	Deionised water	Non venting	White	£25.90
<b>BOT9366</b>	WGW536PML	500	Water	Non venting	White	£26.00

