

Esco Airstream® Class II Biological Safety Cabinets Continued

ESCO
WORLD CLASS. WORLDWIDE.

G-Series

- Tempered glass side wall (UV-absorbing tempered glass, 5mm, colourless and transparent)
- Multi-piece autoclavable work tray
- Dual blower fan system
- Dual exhaust filter

Code	Alt Ref	Model	Nominal Size, ft	Ext Dims, w x d* x h, mm	Int Dims, w x d x h, mm	Price
CAB7522	2010734	AC2-4G8	4	1340 x 753 x 1400	1220 x 580 x 660	£8466.00
CAB7524	2010743	AC2-6G8	6	1950 x 753 x 1400	1830 x 580 x 660	£11817.00

* Without armrest

D-Series

- One piece stainless steel side wall (1.5mm 16 gauge stainless steel, type 304, with 4B finish)
- Single piece stainless steel work tray
- Dual blower fan system
- Dual exhaust filter

Code	Alt Ref	Model	Nominal Size, ft	Ext Dims, w x d* x h, mm	Int Dims, w x d x h, mm	Price
CAB7518	2010733	AC2-4D8	4	1340 x 753 x 1400	1220 x 580 x 660	£8946.00
CAB7520	2010742	AC2-6D8	6	1950 x 753 x 1400	1830 x 580 x 660	£12296.00

* Without armrest



As a manufacturer of laboratory equipment, Esco offers tailored solutions that fit the needs of laboratories in various industries



**SEE OUR
CO₂ INCUBATOR
IN SECTION I**

CelCulture® CO₂ Incubator



SCAN HERE

- SwiftCon™ 90 °C validated moist heat decontamination cycle (15 hours)
- ULPA filtration ISO class 5 air quality and 170L capacity



SCAN HERE



SCAN HERE

FIND OUT MORE



www.escolifesciences.com



@Esco-Great-Britain



ESCO
LIFESCIENCES GROUP

Servicing a Range of Laboratory Equipment by Esco

- Inventory of available parts, technical support, field-based and factory engineers
- Annual Testing/Re-Certification
- Troubleshooting and Parts Replacement
- Repair and Maintenance
- Installation and Operational Qualification