

## Esco Airstream® Plus Class II Microbiological Safety Cabinets Continued



## Accessories

Code	Alt Ref	Description	For use with	Price
<b>CAB7308</b>	5130134	Fixed (28") support stand with levelling feet	4ft cabinets	£618.00
<b>CAB7312</b>	5130171	Fixed (28") support stand with levelling feet	5ft cabinets	£675.00
<b>CAB7316</b>	5130172	Fixed (28") support stand with levelling feet	6ft cabinets	£675.00
<b>CAB7310</b>	5130175	Fixed (34") support stand with levelling feet	4ft cabinets	£618.00
<b>CAB7314</b>	5130176	Fixed (34") support stand with levelling feet	5ft cabinets	£675.00
<b>CAB7318</b>	5130177	Fixed (34") support stand with levelling feet	6ft cabinets	£675.00
<b>CAB1232</b>	5130051	Telescoping (28 - 34") support stand with levelling feet	4ft cabinets	£738.00
<b>CAB1233</b>	5130052	Telescoping (28 - 34") support stand with levelling feet	5ft cabinets	£858.00
<b>CAB1234</b>	5130053	Telescoping (28 - 34") support stand with levelling feet	6ft cabinets	£858.00
<b>CAB1228</b>	5130056	Telescoping (28 - 34") support stand with castors	4ft cabinets	£738.00
<b>CAB1229</b>	5130057	Telescoping (28 - 34") support stand with castors	5ft cabinets	£858.00
<b>CAB1230</b>	5130058	Telescoping (28 - 34") support stand with castors	6ft cabinets	£858.00

## Haier Class II Biological Safety Cabinets

**Haier Biomedical**  
Intelligent Protection of Life Science

Provide an effective and greener solution in operator, product and environmental protection. These biosafety cabinets are BS EN 12469 and TUV certified. Energy efficient Class II microbiological safety cabinet complete with stand, two patented IP44 rated waterproof power sockets, LCD screen, four service ports and UV lamp for sterilisation.

- Energy efficient with DC and EC fans and longer last LED lights which reduce energy consumption, reduce noise and increase air flow uniformity
- Choice of single or two highly efficient HEPA dual out filters reduce the risk of lab pollution
- Maximum comfort - ergonomic design ensures maximum comfort and alleviates operator fatigue
- Certified to BS EN 12469 standard
- Multiple alarm functions including fan turn off alarm after door opening, door open alarm, filter alerts
- Adjustable height stand (680 - 900mm) with hidden mechanism to avoid contamination
- CAB1264 supplied without stand



Code	Alt Ref	Ext Dims, w x d x h, mm	Int Dims, w x d x h, mm	Price
<b>CAB1260</b>	HR900-IIA2	1080 x 845 x 2160	920 x 620 x 650	£8999.00
<b>CAB1262</b>	HR1200-IIA2	1380 x 845 x 2160	1220 x 620 x 650	£9901.00
<b>CAB1264</b>	HR1200-IIA2-D-NS	1380 x 780 x 1460	1310 x 620 x 630	£10481.00
<b>CAB1244</b>	HR1200-IIA2-D	1380 x 780 x 2160	1310 x 620 x 630	£11366.00
<b>CAB1266</b>	HR1500-IIA2	1680 x 845 x 2160	1520 x 620 x 650	£10883.00

## View our Safety &amp; PPE catalogue

To browse our new Safety and PPE catalogue, visit the SLS website:

**[scientificlabs.co.uk/literature](https://scientificlabs.co.uk/literature)** (UK & Northern Ireland)

**[scientificlabs.ie/literature](https://scientificlabs.ie/literature)** (Republic of Ireland)

**[scientificlabs.com/literature](https://scientificlabs.com/literature)** (International)

