

## Rodwell Ambassador Autoclaves



The Ambassador is a heavy-duty autoclave that is capable of meeting the needs of any busy laboratory.

- MP25 controller
- Front loading, cylindrical chamber
- Robust, space-saving design
- Automatic fill, drain and refill
- Scale prevention systems
- Waste discard system
- High speed fan accelerated cooling
- Range of cycle recording systems
- Air removal HEPA filter
- Air ballasting
- Vacuum air removal



Code	Alt Ref	Steam Source	Capacity, L	Chamber Dims, Ø x d, mm	Dims, w x d x h, mm	Power Supply	Price
<b>AUT1301</b>	AMB24	Electronically heated	100	460 x 610	650 x 1160 x 1240	3ph/N 50/60Hz	<b>£15700.00</b>
<b>AUT1297</b>	AMB28	Electronically heated	158	460 x 960	650 x 1310 x 1240	3ph/N 50/60Hz	<b>£16750.00</b>
<b>AUT1298</b>	AMB29	Direct steam	158	460 x 960	650 x 1310 x 1240	1ph/N 50/60Hz	<b>£16750.00</b>

## Options and Accessories

Code	Alt Ref	Description	Price
<b>AUT2007</b>	AC30	Waste discard system	<b>£1244.00</b>
<b>AUT1356</b>	AC2B	Fan assisted cooling	<b>£921.00</b>
<b>AUT2046</b>	AC260B	Compressed air ballasting	<b>£1669.00</b>
<b>AUT1376</b>	AC260C	Oil free compressor for ballasting	<b>£2064.00</b>
<b>AUT2019</b>	AC164	Condense collector	<b>£277.00</b>
<b>AUT1358</b>	AC110	Drain line condenser for plastic drains	<b>£1071.00</b>
<b>AUT2015</b>	AC61	Additional loading basket	<b>£1294.00</b>
<b>AUT1349</b>	AC90	Automatic fill	<b>£1151.00</b>
<b>AUT1351</b>	AC91	Automatic fill/drain for clean steam	<b>£1847.00</b>
<b>AUT1377</b>	SBB1	Under-chamber boiler (c/w autofill)	<b>£3300.00</b>
<b>AUT1379</b>	PHB1	Pre-heated under chamber boiler	<b>£3994.00</b>
<b>AUT1381</b>	AC256	Condensate removal for boiler	<b>£1077.00</b>



**Our service options:**

- Project Design & Management
- Installation
- Qualification
- Repair
- Maintenance
- Calibration
- Validation
- Spare Parts
- Training



## Clarity in Servicing

**Our Knowledge + Your Equipment = Service Tailored to You**

**Email:** [servicecentre@scientific-labs.com](mailto:servicecentre@scientific-labs.com)    **Telephone:** 0115 977 5075