



Esco Ascent™ Ductless Fume Hoods

ESCO
WORLD CLASS. WORLDWIDE.

Ensures superior chemical fume containment and offers a sensible balance of quality and performance features. The green solution – environmentally safe and does not discharge toxic gases to the environment.

- Safe Carbon filtration – complies to all international standards
- Large working area provides user comfort
- Generous carbon filter weight to ensure uniformed filtration
- Elegant and smooth sliding motorised sash window
- Isocide™ powder coating prevents surface contamination
- Grade 304 stainless steel worktop surface that will never rust, chip or generate particles
- Electrogalvanised steel wall for maximum durability
- Lip at the front edge of work surface ensures that any spillage is contained in the work zone
- Sentinel™ microprocessor control system with audible and visual alarms for unsafe conditions
- Optional VOC sensor to monitor chemical concentration within the cabinet work zone
- All cabinet service and filter replacement can be carried out from the front, allowing the cabinet to be placed against wall in the laboratory to save space
- No ductwork or exhaust system required
- VOC sensors, swan neck faucets and electrical outlets are available as factory fitted options – please contact Customer Services for a tailored quotation.
- Secondary carbon/HEPA filter models also available – please enquire

Code	Alt Ref	Model	Size	Int Dims, w x d x h, mm	Ext Dims, w x d x h, mm	Price
CAB7530	2040024	ADC-3B1	3ft/0.9m	965 x 596 x 765	1035 x 736 x 1455	£5611.00
CAB7534	2040042	ADC-4B1	4ft/1.2m	1268 x 596 x 765	1340 x 736 x 1455	£5611.00
CAB7537	2040060	ADC-5B1	5ft/1.5m	1575 x 596 x 765	1645 x 736 x 1455	£6422.00
CAB7550	2040078	ADC-6B1	6ft/1.8m	1880 x 596 x 765	1950 x 736 x 1455	£6153.00

Support Stands with Castor Wheels

- Height: 711mm

Code	Alt Ref	Model	For use with	Price
CAB7544	5130155	SPC-3A0 Gen 2	Fume hoods 1035mm wide	£618.00
CAB7546	5130152	SPC-4A0 Gen 2	Fume hoods 1340mm wide	£618.00
CAB7548	5130167	SPC-5A0 Gen 2	Fume hoods 1645mm wide	£675.00
CAB7552	5130154	SPC-6A0 Gen 2	Fume hoods 1950mm wide	£675.00



FASTER ChemFAST Elite Fume Cupboards

FASTER

ChemFAST belong to the latest generation of ductless fume cupboards using the most innovative molecular filtration technology manufactured by FASTER, in which the choice of materials of construction of the highest quality guarantees a safe working environment together with fume containment for protection from chemicals, vapours and aerosol in the laboratories.

- Microprocessor based monitoring system for a range of data which includes type and code number of the filter being used, installation date, maximum time allowance for filter use and a warning date for its replacement
- Monitoring data: power on/off, light on/off, microprocessor monitoring system checking airflow, pre-filter and filter efficiency, variable speed air regulation, audible and visual alarms alert the operator to low/high airflow, fan failure, filter and pre-filter condition, black-out, gas detector and anemometer failure
- External structure made by cold rolled epoxy painted galvanised for superior resistance against chemical solvents
- Acid resistant PVC rear wall and work surface with special design to drive the liquids spilt in the corners of the surface