

Esco Airstream® Laminar Flow Cabinets Continued



Accessories and Options Continued

Code	Alt Ref	Description	For use with	Price
CAB1232	5130051	Telescoping (28 - 34") support stand with levelling feet	4ft cabinets	£738.00
CAB1233	5130052	Telescoping (28 - 34") support stand with levelling feet	5ft cabinets	£858.00
CAB1234	5130053	Telescoping (28 - 34") support stand with levelling feet	6ft cabinets	£858.00
CAB7304	5130170	Fixed (28") support stand with levelling feet	3ft cabinets	£618.00
CAB7308	5130134	Fixed (28") support stand with levelling feet	4ft cabinets	£618.00
CAB7312	5130171	Fixed (28") support stand with levelling feet	5ft cabinets	£675.00
CAB7316	5130172	Fixed (28") support stand with levelling feet	6ft cabinets	£675.00
CAB7306	5130174	Fixed (34") support stand with levelling feet	3ft cabinets	£618.00
CAB7310	5130175	Fixed (34") support stand with levelling feet	4ft cabinets	£618.00
CAB7314	5130176	Fixed (34") support stand with levelling feet	5ft cabinets	£675.00
CAB7318	5130177	Fixed (34") support stand with levelling feet	6ft cabinets	£675.00

SLS Lab Basics Mini Centrifuge



Ideal for microfiltration and quick spin downs from the walls and caps of microcentrifuge tubes. The compact design saves bench space and features an easy-open clear lid. Simply close the lid for acceleration and open the lid for deceleration, or alternatively use the on/off switch to control.

- Fixed speed
- Easy to use with dual start/stop function
- Compact design
- Whisper quiet operation
- Ideal for quick spin downs of microtubes and PCR tubes
- Universal voltage
- Rotor change without tools
- Accelerates and brakes in seconds
- Supplied with 8 place 1.5/2.0mL microtube rotor and 16 place (2 x 8 strips) PCR rotor



Code	Max Speed, RPM	Max RCF, xg	Dims, w x d x h, mm	Price
SLS1660	7000	2680	160 x 170 x 122	£175.00

SLS Lab Basics my-FUGE Mini Centrifuge



Eliminates the need to change rotors when switching between microtubes and PCR strips. Includes a unique COMBI-Rotor™ allowing the running of 12 x 1.5/2.0mL microtubes and 4 PCR strips simultaneously.

- Fixed speed
- Twice the capacity of most mini centrifuges
- Near silent operation
- Starts and stops with closing/opening of the lid
- Conserves valuable bench space



Code	Max Speed, RPM	Max RCF, xg	Dims, w x d x h, mm	Price
CEN5050	5500	2000	146 x 203 x 108	£538.00

Adapters

Supplied in packs of 6

Code	Description	Price
CEN5054	Adaptor for 0.2mL tubes	£22.00
CEN5052	Adaptor for 0.5mL tubes	£22.00