



Corning® LSE Mini Microcentrifuge

CORNING

A compact, lightweight and easy to use centrifuge. It includes an 8 place rotor that will accept standard 1.5 to 2.0mL microcentrifuge tubes, and adapters are included for 0.5/0.4mL and 0.2mL tubes. Also included is a 4 place rotor for PCR strip tubes.

- Fixed speed
- Close the lid to start, press the lid release button to stop
- Quick release rotor system eliminates the need for tools to exchange rotors
- Electronic braking provides rapid deceleration for more efficient sample prep
- Compact, rugged design withstands heavy use in crowded lab environment
- Dual lid switch provides safe operation
- Ergonomic design allows for comfortable handling

Code	Alt Ref	Max Speed, RPM	Max RCF, xg	Dims, w x d x h, mm	Price
CEN1016	6770	6000	2000	148 x 174 x 122	£216.00



Heathrow Scientific Sprout® Plus Mini Centrifuge

 Heathrow Scientific
design | service | innovation

With an ergonomic design and safety measures, the Sprout® Plus is easy to use and brings increased efficiency to the lab. It is ideal for microfilter cell separations and HPLC samples, with adapters to accommodate all major microtube and PCR tube styles (1.5/2.0, 0.5, 0.2mL) for customised use.

- Tool free SnapSpin rotors and quick spin options
- Locking lid designed to fit in the palm of the hand for easy repeat opening and quick tube extraction
- Large lid size accommodates longer and skirted tubes and tall tube lids
- Braking function quickly stops spinning when powered off or when lid is opened
- Weighted rubberised base eliminates potential movement and dampens vibration and noise
- Maximum speed: 6000RPM
- Maximum RCF: 2000xg
- Dims, L x w x h: 153 x 128 x 109mm
- Supplied with 6 place 1.5/2.0mL tube rotor, 16 place 0.2mL tube rotor (single or 2 x 8 strip), 6 x 0.2mL tube adapters, 6 x 0.5mL tube adapters, 4 x lid o-rings, storage case for rotors and adapters, tube rack, and power adapter with 4 interchangeable plugs

Code	Alt Ref	Colour	Price
CEN2114	120610	Green	£331.00
CEN2116	120611	Blue	£331.00



Benchmark StripSpin™ Mini Centrifuge

 Benchmark
Scientific

Designed for centrifugation of 12 position PCR strips. With a fixed speed, this centrifuge is perfect for quick spin downs. Simply close the lid and the unit quickly ramps up. Open the lid, and the rotor quickly decelerates for easy removal of samples.

- Compact footprint
- Near silent operation
- Starts and stops with closing/opening of the lid
- Supplied with unique rotor that holds 12 position PCR strips
- Optional rotor available for 8 x 1.5/2.0mL tubes

Code	Alt Ref	Max Speed, RPM	Max RCF, xg	Dims, w x d x h, mm	Price
CEN1714	C1248-E	5500	2000	140 x 200 x 112	£553.00

Accessories

Code	Alt Ref	Description	Price
CEN1715	C1008-ROT	Optional rotor for 8 x 1.5/2.0mL tubes	£53.00

