



Eppendorf 5425 Centrifuges

eppendorf

24 place microcentrifuges engineered to be extremely quiet. Short runs up to 12 minutes are whisper quiet at an average noise level of 45dB(A) and maximum speed of up to 21300xg. This is sufficient for centrifugation steps of most routinely used DNA/RNA purification kits. During longer runs the cooling fan works at full capacity to dissipate the heat for the remainder of the run.

- Aerosol tight Eppendorf QuickLock® lids for simple and safe handling
- Outstanding software features (e.g. adjustable acceleration and braking ramps, 3 program buttons, short function without continuous pressing of the key)
- OptiBowl® design for quiet operation, even without a rotor lid
- Soft touch one finger lid closure for ergonomic operation
- Timer: 10s to 9:59hr, with continuous run function
- Dims, w x d x h: 240 x 390 x 240mm
- Speed range: 100 – 15060RPM
- Max RCF: 21300xg (with fixed angle rotor); 3217xg (with swing bucket rotor)

Code	Description	Supplied with	Capacity	Price
E5405000360	5425 centrifuge with keypad	—	—	£1797.00
E5405000166	5425 centrifuge with keypad with FA-24x2 rotor	FA-24x2 rotor	24 x 1.5/2.0mL tubes	£2234.00
E5405000760	5425 centrifuge with rotary knobs	—	—	£1797.00
E5405000565	5425 centrifuge with rotary knobs	FA-24x2 rotor	24 x 1.5/2.0mL tubes	£2234.00

Rotors and Adapters

Code	Description	Capacity	Max Speed, RPM	Max RCF, xg	Pack	Price
E549550008	FA-10x5 rotor with aerosol tight QuickLock® lid	10 x 5.0mL tubes (snap and screw cap)	15060	21300	Each	£738.00
E5495500006	FA-24x2 rotor with aerosol tight QuickLock® lid	24 x 1.5/2.0mL tubes	15060	21300	Each	£579.00
E5495503005	FA-24x2-PTFE rotor with aerosol tight QuickLock® lid	24 x 1.5/2.0mL tubes	15060	21300	Each	£738.00
E5495508007	FA-18x2-Kit rotor with aerosol tight QuickLock® lid	18 x 1.5/2.0mL tubes or spin columns	15060	18565	Each	£460.00
E5495510001	F-32x0.2-PCR rotor with lid	32 x 0.2mL PCR tubes or 4 x 8 PCR strips	15060	18134	Each	£422.00
E5495512004	S-96x0.2-PCR rotor with buckets	96 x 0.2mL PCR tubes or 12 x 8 PCR strips	6000	3217	Each	£579.00
E5425715005	Adapter for all 1.5/2.0mL rotors	1 x 0.2mL PCR tube (max Ø 6mm)	—	—	6	£12.95
E5425716001	Adapter for all 1.5/2.0mL rotors	1 x 0.5mL tube or BD Microtainer® (max Ø 8mm)	—	—	6	£12.95
E5425717008	Adapter for all 1.5/2.0mL rotors	1 x 0.4mL tube (max Ø 6mm)	—	—	6	£12.95
E5820768002	Adapter for all 5.0mL rotors	1 x 1.5 - 2.0mL vessel	—	—	4	£31.00
E5820769009	Adapter for all 5.0mL rotors	1 x cryo tube	—	—	4	£30.60
E5820770007	Adapter for all 5.0mL rotors	1 x HPLC tube	—	—	4	£30.60

A world of products in one place...

- Over 500,000 products
- Easy ordering capability
- Detailed product information
- Promotions and offers

