

Eppendorf 5702/5702 R/5702 RH Centrifuges Continued

eppendorf

Code	Description	Capacity	Max RCF, xg	Max Speed, RPM	Pack	Price
E5702721000	Aerosol tight caps 100mL round bucket in A-4-38 rotor	—	—	—	2	£54.00
E5702717002	Adapters for 100mL round bucket in A-4-38 rotor	25mL conical tubes	—	—	2	£83.00
E5702718009	Adapters for 100mL round bucket in A-4-38 rotor	1 x 85/100mL round bottom tubes	—	—	2	£83.00
E5702724009	Adapters for 100mL round bucket in A-4-38 rotor	4 x 9 - 15mL round bottom tubes	—	—	2	£149.00
E5702745006	Adapters for 100mL round bucket in A-4-38 rotor	4 x 1.5 - 2mL reaction vessels	—	—	2	£154.00
E5702746002	F-45-24-11 rotor	24 x 1.5/2.0mL tubes	1770	4400	Each	£354.00
E5702747009	F-45-18-17 cryo rotor	1 - 2mL cryo container (max 13mm Ø, use with adapter) and tubes with screw cap (max 17mm Ø, no adapter required)	1770	4400	Each	£406.00

Eppendorf 5910 Ri Refrigerated Centrifuges

eppendorf

Combine high capacity and performance in a compact and ergonomic footprint. The large 7 inch VisioNize® touch interface with an intuitive user experience provides fast and error free operation, ensuring efficient and reproducible runs. Additionally, the documentation function and unique user management support your quality control with multiple authorisation levels, password protection, and documentation of all runs and user activities.

- User management with multiple authorization levels and documentation features supporting GxP requirements
- High centrifugation speed of up to 22132xg (14000RPM)
- Aerosol tight Eppendorf QuickLock® caps and lids for fast and ergonomic lid locking
- Quiet operation to improve your work environment
- Centrifuge lid with soft touch lid closure
- Temperature range from -11 to 40°C
- Automatic rotor recognition and imbalance detection for maximum operational safety



Code	Description	Supplied with	Price
E5943000061	5910 Ri refrigerated centrifuge without rotor	—	£8582.00
E5943000157	5910 Ri refrigerated centrifuge	Rotor FA-6x50	£9944.00
E5943000181	5910 Ri refrigerated centrifuge	Rotor S-4xUniversal with universal buckets 1000mL bottles and 1000mL adapters	£11428.00
E5943000266	5910 Ri refrigerated centrifuge	Rotor S-4xUniversal with universal buckets adapters for 5/15/50mL and 175-250mL bottles and plates	£11429.00
E5943000361	5910 Ri refrigerated centrifuge	Rotor S-4xUniversal with universal buckets and adapters for 5/15/50mL/plates	£11340.00
E5943000561	5910 Ri refrigerated centrifuge	Rotor S-4x750 with round buckets and adapters for 5/15/50mL/plates	£10824.00
E5943000912	5910 Ri refrigerated centrifuge	Rotor S-4xUniversal with universal buckets adapters for 5/15/50mL and 175-250mL bottles and plates and 1000mL bottles 1000mL adapters	£12075.00

Rotors and Accessories

Code	Description	Capacity	Max RCF, xg	Max Speed, RPM	Pack	Price
E5895201008	S-4xUniversal rotor without buckets	—	4347	4500	Each	£1312.00
E5895200001	S-4xUniversal rotor with 4 x universal buckets	—	4347	4500	Each	£2946.00
E5895203000	Universal buckets for S-4xuniversal rotor	—	—	—	4	£1663.00
E5895202004	Universal buckets for S-4xuniversal rotor	—	—	—	2	£854.00
E5895181007	S-4x400 rotor without buckets	—	5263	5100	Each	£1044.00
E5895180000	S-4x400 rotor with 4 x 400mL round buckets	—	5263	5100	Each	£2186.00
E5895182003	Round buckets for S-4x400 rotor	400mL	—	—	4	£1331.00
E5895183000	Round buckets for S-4x400 rotor	400mL	—	—	2	£718.00
E5910701004	Adapter for S-4x400 rotor	4 x 50mL conical tubes	—	—	2	£207.00