



## Eppendorf Cell Culture and Clinical Centrifuge Packages

eppendorf

The Eppendorf centrifuge packages offer a selection of centrifuges supplied with rotors and adapters to meet specific needs.

**Cell Culture Packages:** optimised for swing bucket rotor applications in standard 15mL and 50mL Falcon® tubes and the new swing bucket rotors increase the capacity of Centrifuges 5804/5804 R to 4 x 250mL and Centrifuges 5810/5810 R to 4 x 750mL. They include:

- Centrifuge 5702, 5804/ 5804 R or 5810/ 5810 R
- Swing bucket rotor including buckets
- 4 x adapters for 15mL Falcon® tubes (2 sets of 2)
- 4 x adapters for 50mL Falcon® tubes (2 sets of 2)

**Clinical Packages:** optimised for Centrifuge 5702 R with adapters for 13mm and 16mm blood collection tubes or urine tubes. They include:

- Centrifuge 5702 R
- Swing bucket rotor A-4-38 including round buckets
- 4 x adapters for 13mm blood collection tubes or urine tubes (2 sets of 2)
- 4 x adapters for 16mm blood collection tubes or urine tubes (2 sets of 2)

| Cell Culture Centrifuge Packages | Number of Vessels per Adapter/Rotor 15mL Conical | Number of Vessels per Adapter/Rotor 50mL Conical |
|----------------------------------|--------------------------------------------------|--------------------------------------------------|
| 5702/5702R                       | 2/8                                              | 1/4                                              |
| 5804/5804R + rotor A-4-44        | 4/16                                             | 1/4                                              |
| 5804/5804R + rotor S-4-72        | 8/32                                             | 4/16                                             |
| 5810/5810R + rotor A-4-62        | 9/36                                             | 4/16                                             |
| 5810/5810R + rotor A-4-81        | 12/48                                            | 5/20                                             |
| 5810/5810R + rotor S-4-104       | 14/56                                            | 7/28                                             |
| 5810/5810R + rotor S-4-104       | 23/92 13mm tubes                                 | 20/80 16mm tubes                                 |

| Code        | Description                    | Supplied with                                                                                        | Price    |
|-------------|--------------------------------|------------------------------------------------------------------------------------------------------|----------|
| E5702000365 | 5702 centrifuge                | A-4-38 swing-bucket rotor, round buckets, 15mL/50mL conical adapters                                 | £2579.00 |
| E5703000360 | 5702 R refrigerated centrifuge | A-4-38 swing-bucket rotor, round buckets, adapters for 13mm and 16mm blood collection or urine tubes | £4629.00 |
| E5804000360 | 5804 centrifuge                | A-4-44 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters                           | £4420.00 |
| E5804000565 | 5804 centrifuge                | S-4-72 rotor with adapters for 15mL/50mL conical tubes                                               | £5312.00 |
| E5805000365 | 5804 R refrigerated centrifuge | A-4-44 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters                           | £6734.00 |
| E5805000665 | 5804 R refrigerated centrifuge | S-4-72 rotor, adapters for 15mL/50mL conical tubes                                                   | £7400.00 |
| E5810000365 | 5810 centrifuge                | A-4-62 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters                           | £5790.00 |
| E5810000460 | 5810 centrifuge                | A-4-81 swing-bucket rotor, rectangular buckets, 15mL/50mL conical                                    | £6510.00 |
| E5810000560 | 5810 centrifuge                | S-4-104 rotor, adapters for 15mL/50mL conical tubes                                                  | £6972.00 |
| E5811000360 | 5810 R refrigerated centrifuge | A-4-62 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters                           | £8433.00 |
| E5811000465 | 5810 R refrigerated centrifuge | A-4-81 swing-bucket rotor, rectangular buckets, 15mL/50mL conical adapters                           | £9152.00 |
| E5811000660 | 5810 R refrigerated centrifuge | S-4-04 rotor, adapters for 15mL/50mL conical tubes                                                   | £9614.00 |
| E5810000665 | 5810 centrifuge                | S-4-104 swing-bucket rotor, adapters for 13-16mm round-bottom tubes                                  | £6971.00 |
| E5811000865 | 5810 R refrigerated centrifuge | S-4-104 swing-bucket rotor, adapters for 13/16mm round bottom tubes                                  | £9614.00 |