

## Thermo Scientific™ SL 1 Plus Series Centrifuges Continued

**Thermo**  
SCIENTIFIC

## Rotors and Accessories

Code	Description	Capacity	Max Speed, RPM	Max RCF, xg	Pack	Price
<b>75003181</b>	TX-400 swinging bucket rotor (without buckets)	4 x 400mL	5000	4696	Each	<b>£1561.00</b>
<b>75003655</b>	Round buckets 400mL for TX-400 rotor	—	—	—	4	<b>£1398.00</b>
<b>75003656</b>	ClickSeal™ biocontainment lids for TX-400 rotor buckets	—	—	—	4	<b>£387.00</b>
<b>75003680</b>	Adapters for TX-400 rotor	19 x 5/7mL vacutainers	—	—	4	<b>£393.00</b>
<b>75003681</b>	Adapters for TX-400 rotor	14 x 10/15mL tubes	—	—	4	<b>£393.00</b>
<b>75003682</b>	Adapters for TX-400 rotor	9 x 15mL conical	—	—	4	<b>£393.00</b>
<b>75003683</b>	Adapters for TX-400 rotor	4 x 50mL conical	—	—	4	<b>£393.00</b>
<b>75003706</b>	Adapters for TX-400 rotor	3 x 25mL universal	—	—	4	<b>£417.00</b>
<b>75003788</b>	Adapters for TX-400 rotor	1 x 250mL oak ridge or 200mL/175mL conical bottles	—	—	4	<b>£394.00</b>
<b>75003300</b>	H-Flex swinging bucket rotor (without buckets)	16 x 50mL	4700	4297	Each	<b>£2477.00</b>
<b>75003301</b>	Rectangular buckets for H-Flex rotor	—	—	—	2	<b>£1566.00</b>
<b>75003302</b>	Biocontainment lids for H-Flex rotor bucket	—	—	—	2	<b>£353.00</b>
<b>75003303</b>	Adaptors for H-Flex rotor	4 x 50mL conical tubes	—	—	2	<b>£718.00</b>
<b>75003304</b>	Adaptors for H-Flex rotor	9 x 15mL conical tubes	—	—	2	<b>£718.00</b>
<b>75003305</b>	Adaptors for H-Flex rotor	33 x 10/12mL tubes	—	—	2	<b>£718.00</b>
<b>75003306</b>	Adaptors for H-Flex rotor	42 x 5/7mL tubes	—	—	2	<b>£718.00</b>
<b>75003307</b>	Adaptors for H-Flex rotor	5 x microplates or 2 x deep well plates	—	—	2	<b>£718.00</b>
<b>75003698</b>	Fiberlite™ F15-6 x 100y fixed angle rotor	6 x 100mL	15000	24652	Each	<b>£2629.00</b>
<b>75003103</b>	Adapter for F15-6 x 100y rotor	6 x 50mL conical tubes	—	—	Each	<b>£75.00</b>
<b>75003624</b>	M-20 swinging bucket microplate rotor (without carriers)	2 x 3 standard microplates, 2 x 1 deepwell	4000	2272	Each	<b>£1714.00</b>
<b>75003625</b>	Sealed microplate carriers (tray and lids) for M-20 rotor	—	—	—	2	<b>£970.00</b>
<b>75002011</b>	ClickSeal™ replacement lids for M-20 rotor	—	—	—	4	<b>£193.00</b>



2 years warranty, 3 years extra on drive and refrigerant

## Thermo Scientific™ X4/X4R Multifuge Centrifuges

**Thermo**  
SCIENTIFIC

Thermo Scientific™ centrifuges offering continuous innovation with bio containment, design, performance, ergonomics & value. The X4/X4R series is suitable for all sample preparation specifically cell culture applications, clinical protocols and microplate processing and other separation needs.

- Touch screen user interface with lighthouse mode
- Auto lock rotor technology
- Scheduling features
- Bio containment click seal lids
- Carbon fibre rotors, improved lid design and reach in
- Step programming
- Max speed 15200RPM (25000xg) with fixed angle rotor

Code	Description	Price
<b>75009900</b>	Multifuge X4 Pro	<b>£6822.00</b>
<b>75009515</b>	Multifuge X4R Pro	<b>£11063.00</b>
<b>75016067</b>	Multifuge X4 Pro cell culture package	<b>£11645.00</b>
<b>75016066</b>	Multifuge X4R Pro cell culture package	<b>£15686.00</b>