

Thermo Scientific™ SL 1 Plus Series Centrifuges Continued

Rotors and Accessories

Code	Description	Capacity	Max Speed, RPM	Max RCF, xg	Pack	Price
75003181	TX-400 swinging bucket rotor (without buckets)	4 x 400mL	5000	4696	Each	£1561.00
75003655	Round buckets 400mL for TX-400 rotor	—	—	—	4	£1398.00
75003656	ClickSeal™ biocontainment lids for TX-400 rotor buckets	—	—	—	4	£387.00
75003680	Adapters for TX-400 rotor	19 x 5/7mL vacutainers	—	—	4	£393.00
75003681	Adapters for TX-400 rotor	14 x 10/15mL tubes	—	—	4	£393.00
75003682	Adapters for TX-400 rotor	9 x 15mL conical	—	—	4	£393.00
75003683	Adapters for TX-400 rotor	4 x 50mL conical	—	—	4	£393.00
75003706	Adapters for TX-400 rotor	3 x 25mL universal	—	—	4	£417.00
75003788	Adapters for TX-400 rotor	1 x 250mL oak ridge or 200mL/175mL conical bottles	—	—	4	£394.00
75003300	H-Flex swinging bucket rotor (without buckets)	16 x 50mL	4700	4297	Each	£2477.00
75003301	Rectangular buckets for H-Flex rotor	—	—	—	2	£1566.00
75003302	Biocontainment lids for H-Flex rotor bucket	—	—	—	2	£353.00
75003303	Adaptors for H-Flex rotor	4 x 50mL conical tubes	—	—	2	£718.00
75003304	Adaptors for H-Flex rotor	9 x 15mL conical tubes	—	—	2	£718.00
75003305	Adaptors for H-Flex rotor	33 x 10/12mL tubes	—	—	2	£718.00
75003306	Adaptors for H-Flex rotor	42 x 5/7mL tubes	—	—	2	£718.00
75003307	Adaptors for H-Flex rotor	5 x microplates or 2 x deep well plates	—	—	2	£718.00
75003698	Fiberlite™ F15-6 x 100y fixed angle rotor	6 x 100mL	15000	24652	Each	£2629.00
75003103	Adapter for F15-6 x 100y rotor	6 x 50mL conical tubes	—	—	Each	£75.00
75003624	M-20 swinging bucket microplate rotor (without carriers)	2 x 3 standard microplates, 2 x 1 deepwell	4000	2272	Each	£1714.00
75003625	Sealed microplate carriers (tray and lids) for M-20 rotor	—	—	—	2	£970.00
75002011	ClickSeal™ replacement lids for M-20 rotor	—	—	—	4	£193.00



2 YEARS



2 years warranty, 3 years extra on drive and refrigerator

Thermo Scientific™ X4/X4R Multifuge Centrifuges

Thermo Scientific™ centrifuges offering continuous innovation with bio containment, design, performance, ergonomics & value. The X4/X4R series is suitable for all sample preparation specifically cell culture applications, clinical protocols and microplate processing and other separation needs.

- Touch screen user interface with lighthouse mode
- Auto lock rotor technology
- Scheduling features
- Bio containment click seal lids
- Carbon fibre rotors, improved lid design and reach in
- Step programming
- Max speed 15200RPM (25000xg) with fixed angle rotor

Code	Description	Price
75009900	Multifuge X4 Pro	£6822.00
75009515	Multifuge X4R Pro	£11063.00
75016067	Multifuge X4 Pro cell culture package	£11645.00
75016066	Multifuge X4R Pro cell culture package	£15686.00