

## Hettich EBA 200/200 S Compact Centrifuges Continued



Code	Alt Ref	Description	Price
<b>CEN0171</b>	GB-C1800	EBA 200 centrifuge with 8 x 15mL angle rotor	£1405.00
<b>CEN0172</b>	GB-C1802	EBA 200 S refrigerated centrifuge with 8 x 15mL angle rotor	£2609.00

## Accessories

Code	Alt Ref	Description	Pack	Price
<b>CEN0174</b>	GB-A6305	Adapter for 8 x 4mL tubes	Each	£34.10
<b>CEN0175</b>	GB-A1054-A	Adapter for 8 x 5mL tubes	Each	£21.00
<b>CEN0176</b>	GB-A1058	Adapter for 8 x 4 - 7mL urine tubes	Each	£35.00
<b>CEN5040</b>	GB-A1063-8	Adapter for 1 x 1.5 - 2.0mL tubes	8	£97.00

## Hettich ROTANTA 460/460 R Centrifuges



The ROTANTA 460 family of centrifuges is an extraordinary line that meets a wide range of demands. They feature high capacity, versatility and performance. Their state of the art technology and high safety standards ensure smooth running for any laboratory application.

- Brushless induction motor
- Medical Device according to directive 93/42/EC
- 28 individual acceleration and deceleration stages
- Easy operation with keypad and control knob
- 98 program memories for more individuality
- Metal housing and lid with viewing port, powered one hand lid lock, dropping protection and emergency lock release
- Automatic rotor recognition
- Imbalance monitoring and switch off



2  
YEARS



	ROTANTA 460	ROTANTA 460 R
Cooling	No	Yes
Temperature Control Range, °C	—	-20 to +40
Max Speed, RPM	15000	15000
Max RCF, g	24400	24400
Dims, w x d x h, mm	554 x 706 x 456	770 x 706 x 456
Weight, kg	101	141

Code	Alt Ref	Description	Price
<b>CEN0701</b>	GB-C5650	ROTANTA 460 centrifuge (without rotor)	£7632.00
<b>CEN0703</b>	GB-C5660	ROTANTA 460 R refrigerated centrifuge (without rotor)	£11086.00

## Rotors and Accessories

Code	Alt Ref	Description	Capacity	Max Speed, RPM	Max RCF, xg	Price
<b>CEN0752</b>	GB-A5645	Angle rotor, 6 place, 25°	6 x 250mL tubes	8500*	11228**	£3176.00
<b>CEN0754</b>	GB-A5615	Angle rotor, 6 place, 35°	6 x 94mL tubes	11500	18038	£2936.00
<b>CEN0756</b>	GB-A4489-A	Angle rotor, 30 place, 45°	30 x 0.2-2mL tubes	15000	24400	£1140.00
<b>GBA5699</b>	GB-A5699-R	Swing-out rotor, 4 place, 90°	—	5063	4600	£1741.00
<b>GBA4880</b>	GB-A4880	Bucket for GB-A5699-R rotor	—	—	—	£401.00
<b>GBA4883</b>	GB-A4883	Bioseal lid for GB-A5699-R rotor bucket	—	—	—	£76.00
<b>A4832</b>	GB-A4832	Adaptor for GB-A5699-R rotor	19 x 1-7mL blood tubes	—	—	£155.00
<b>CEN0660</b>	GB-A4833	Adaptor for GB-A5699-R rotor	19 x 4-15mL tubes	—	—	£155.00
<b>CEN0666</b>	GB-A4839	Adaptor for GB-A5699-R rotor	14 x 12/15mL conical (Falcon type) tubes	—	—	£155.00
<b>A5647</b>	GB-A5647	Adaptor for GB-A5699-R rotor	7 x 50mL conical (Falcon type) tubes	—	—	£170.00
<b>CEN0662</b>	GB-A4834	Adaptor for GB-A5699-R rotor	7 x 25/30mL conical (Falcon type) tubes	—	—	£246.00