

Whatman™ 4 Chr Chromatography Paper

cytiva
Whatman

- Ideal for routine and/or repetitive chromatography when loadings are relatively low
- Manufactured and tested specifically for chromatographic techniques to ensure the wicking capability and uniformity of capillary action that are important in chemical separations
- Very suitable where speed is important and very high resolution is not required
- Pure cellulose produced entirely from the highest quality cotton linters with no additives of any kind
- Also widely used in protein and nucleic acid blotting
- Smooth surface with 0.21mm thickness
- Flow rate 180mm/30min
- Fastest of the thin papers



Code	Alt Ref	Dims, w x l, mm	Pack	Price
CHR1150	3004-917	460 x 570	100	£219.00

SLS Select HPLC Grade Water




Code	Pack	Price
CHE3882	2.5L	£30.00

CHROMATOGRAPHY PORTFOLIO



**SCIENTIFIC
LABORATORY
SUPPLIES**

SLS are partnered with some of the most respected manufacturers to provide a full workflow solution for all your chromatography requirements.



Accucore™
PepMap™
TraceGOLD™
ProPac™
DNAPac™
SMART Digest™
Ascentis® Express
SOLA™
SureSTART™
Superspher®
BioBasic™
Astec®
Sequant®
Hypercarb™
EASY-Spray™
Discovery®
TSKgel®
Ascentis®
Omnifit®
Acclaim™
Purosphere™
HyperSep™
apHera™
GlycanPac™
Titan™ C18
Telos®
Discovery® BIO
SUPELCOGEL™
SeQuant® ZIC-HILIC™
LiChrosorb®
TOYOPEARL®
LiChrospher®
Bioshell™
MabPac™
QuEChERS
LinerGOLD™
Synchronis™
Hypertime
SUPELCOSIL™

Please contact SLS if you are looking for any of these market leading brands