



Tristel Enzystel Triple Enzyme Instrument Cleaner

Tristel

- For the manual cleaning of endoscopes and surgical instruments including soaking in ultrasonic baths
- Effectively dissolves and removes proteinaceous and organic matter, such as blood, saliva and mucus, from the surfaces of instruments and equipment

| Code | Alt Ref | Description | Pack | Price |
|----------------|---------|---------------------------------|--------|----------------|
| TRI2010 | MD315 | Concentrate, self dosing bottle | 6 x 1L | £192.00 |



Tristel Medistel High Level Instrument Disinfectant

Tristel

- For immersible instruments and equipment
- Based on quaternary ammonium compounds
- Primarily designed for open sink use but is also suitable for use in ultrasonic baths
- Effective against bacteria, fungi, yeast and viruses

| Code | Alt Ref | Description | Pack | Price |
|----------------|---------|---------------------------------|--------|----------------|
| TRI1375 | MD311 | Concentrate, self dosing bottle | 6 x 1L | £220.00 |



Tristel Airstel High Level Airborne Surface Disinfectant

Tristel

- Single shot, ready to use airborne disinfectant
- Designed for the bactericidal, fungicidal and virucidal disinfection of inaccessible surfaces
- Unique formulation incorporates a natural cinnamon oil extract which is especially effective against *Aspergillus brasiliensis*
- Designed for terminal disinfection following contact surface disinfection with Tristel Fuse for Surfaces
- Easy to use, with no assembly
- Effective against microorganisms of concern, including coronaviruses

| Code | Alt Ref | Pack | Price |
|----------------|---------|-------|---------------|
| TRI2500 | AG001 | 50mL | £18.80 |
| TRI2502 | AG003 | 300mL | £37.10 |

Can't find what you're looking for?

Let us know! Our team are available to help. Contact your local Territory Sales Manager to discuss your requirements.



**SCIENTIFIC
LABORATORY
SUPPLIES**

