

Tristel SHOT

Tristel

- Designed for the sporidical disinfection of environmental surfaces and non-porous materials in laboratories
- One SHOT generates 500mL of chlorine dioxide working solution
- Effective against bacterial spores, mycobacteria, viruses, fungi, yeast and vegetative bacteria in a contact time of 5 minutes
- Each capsule contains concentrates of base (citric acid) and activator (sodium chlorite), which are released into the SHOT bottle after it has been pre-filled with water
- Small, light and innovative, designed with the environment in mind
- One SHOT delivers 100 x 5mL doses of working solution which equates to 100 surface wipes, all dosed from a system of reusable and recyclable bottles

Code	Description	Pack	Price
CLE2024	Capsule	40 capsules	£431.00
CLE2026	500mL bottle	Each	£3.05
CLE2028	5mL dosing head	Each	£1.68



Rely+On™ Virkon™ Disinfectant

- Unique broad spectrum formulation - no other disinfectant has the same powerful composition or extensive portfolio of performance and safety testing data
- Combines application flexibility with broad spectrum efficacy, on hard surfaces, and in the face of organic challenge
- Independently proven highly effective against over 100 strains of virus in 22 viral families, over 400 strains of bacteria, and over 60 strains of fungi and yeast using a wide variety of contact times, temperatures and organic challenge levels
- Unrivalled spectrum
- Built in colour indicator of strength
- Resistant to inactivation by organic material
- Contains a surfactant - combines cleaning and disinfection in one operation
- Avoids the problems of pre-cleaning which spreads infection
- Tablets are convenient to store and easy to handle

Code	Type	Pack	Makes	Price
CLE1550	Powder	50 x 50g sachets	1 sachet = 5L disinfectant	£210.00
CLE1552	Powder	6 x 500g	1 bottle = 50L disinfectant	£204.00
CLE1554	Powder	5kg	1 drum = 500L disinfectant	£180.00
33012	Tablets	20 x 10 tablets	10 x 5g tablets = 5L disinfectant	£120.00
CLE1556	Tablets	12 x 50 tablets	50 x 5g tablets = 25L disinfectant	£360.00



Virkon™ S Disinfectant

- Developed to specifically address the practical biosecurity requirements of modern day farming and livestock production
- Powerful capabilities of peroxygen based chemistry in the formulation provide a broad spectrum of efficacy with specific focus against viral disease causing organisms
- Effective on porous surfaces, in hard water, at low temperatures and in the presence of organic challenge
- Proven to kill over 500 strains of viruses, bacteria, and fungi and yeasts
- Effective against Newcastle disease virus, avian influenza virus, *Salmonella*, and *Campylobacter*
- Powerful, fast acting, flexible, multipurpose biosecurity disinfectant
- Gold standard boot dip disinfectant

Code	Type	Pack	Makes	Price
CLE1180	Powder	50 x 50g sachets	1 sachet = 5L disinfectant	£116.00
CLE1194	Tablets	12 x 50 tablets	50 x 5g tablets = 25L disinfectant	£165.00

