



Klercide 70/30 IPA



- 70% IPA blended with deionised water
- Alcohol blend is 0.2µm filtered then filled and double bagged in a Grade C cleanroom before being terminally sterilised by gamma irradiation
- Supplied sterile
- Available as either spray system fitted with the SteriShield Delivery System (SDS) which ensures the sterility of the contents is preserved throughout use or as impregnated pouch wipes
- Consistency across the range with the wipe containing the same formulation as the fluid in use
- Irradiated at no less than 25kGy to ensure spore free alcohol
- SDS gives an in use shelf life of 3 months
- Pouch wipes are low particulate non-woven material, supplied in a convenient resealable pouch
- Pouch wipes feature Klartear to enable ease of opening even when wearing gloves

Code	Description	Pack	Price
CLE1026	Spray	6 x 1L	£188.00
CLE1028	Wipes	20 x 15 wipes	£384.00



Klercide Low Residue Quat



- Very low residue quaternary ammonium compound and Water for Injection (WFI).
- Biocide blend is 0.2µm filtered and sterile filled under Grade A laminar flow in a Grade B cleanroom
- Ready to use daily disinfectant with low residue accumulation effect
- Supplied sterile
- Available as either spray system fitted with the SteriShield Delivery System (SDS) which ensures the sterility of the contents is preserved throughout use or as impregnated pouch wipes
- Consistency across the range with the wipe containing the same formulation as the fluid in use
- SDS gives an in use shelf life of 3 months
- Pouch wipes are sterile, heavy weight pre-laundered ultrasonically sealed edge 100% polyester wipes, supplied in a convenient resealable pouch
- Pouch wipes feature Klartear to enable ease of opening even when wearing gloves

Code	Description	Pack	Price
CLE1066	Spray	6 x 1L	£580.00
CLE1068	Wipes	20 x 10 wipes	£571.00