

## Klercide Neutral Detergent

**ECOLAB®**

- Neutral detergent is 0.2µm filtered then filled and double bagged in a Grade C cleanroom before being terminally sterilised by gamma irradiation.
- Ready to use detergent for small and large surfaces
- Supplied sterile
- Available as either spray system fitted with the SteriShield Delivery System (SDS) which ensures the sterility of the contents is preserved throughout use or as impregnated pouch wipes
- Consistency across the range with the wipe containing the same formulation as the fluid in use
- Irradiated at no less than 25kGy to ensure spore free alcohol
- SDS gives an in use shelf life of 3 months
- Pouch wipes are low particulate non-woven material, supplied in a convenient resealable pouch
- Pouch wipes feature Klertear to enable ease of opening even when wearing gloves

Code	Description	Pack	Price
<b>CLE1070</b>	Spray	6 x 1L	<b>£343.00</b>
<b>CLE1072</b>	Wipes	20 x 15 wipes	<b>£492.00</b>



## Klercide Sporicidal Active Chlorine

**ECOLAB®**

- Terminally sterilised, by gamma irradiation, active chlorine blend
- Biocide blend is 0.2µm filtered and double bagged in a Grade C cleanroom
- Ready to use highly effective rotational sporicide for small and large surfaces
- Supplied sterile
- Available as a trigger spray
- Has excellent sporicidal efficacy

Code	Pack	Price
<b>CLE1074</b>	6 x 946mL	<b>£343.00</b>



## Klercide Sporicidal Enhanced Peroxide

**ECOLAB®**

- Enhanced blend of 1.5% hydrogen peroxide
- Biocide blend is 0.2µm filtered then sterile filled and double bagged under Grade A laminar flow.
- Ready to use, safer rotational sporicide for small and large surfaces
- Supplied sterile
- Available as either spray system fitted with the SteriShield Delivery System (SDS) which ensures the sterility of the contents is preserved throughout use or as impregnated pouch wipes
- Consistency across the range with the wipe containing the same formulation as the fluid in use
- SDS gives an in use shelf life of 3 months
- Pouch wipes are low particulate non-woven material, supplied in a convenient resealable pouch
- Pouch wipes feature Klertear to enable ease of opening even when wearing gloves

Code	Description	Pack	Price
<b>CLE1076</b>	Spray	6 x 1L	<b>£486.00</b>
<b>CLE1078</b>	Wipes	20 x 15 wipes	<b>£622.00</b>

