

## Simport® Scientific SecuriTainer™ I Containers



- Designed for collection, transport and storage of specimens
- Shatter resistant which eliminates issues of leakage and evaporation
- Stackable for easy, safe storage and
- Manufactured from virgin polypropylene.
- Magenta polyethylene lids are ribbed for easy opening
- Almost transparent to allow specimens to be viewed without opening.
- Unique tamper evident screw cap using a tab locking mechanism ensures your peace of mind during transport or storage situations where someone might have manipulated the specimen without your prior knowledge
- 95kpa compliant
- Non sterile
- Functional temperature range: -90 to 100°C



| Code           | Alt Ref    | Capacity, mL | Pack | Price   |
|----------------|------------|--------------|------|---------|
| <b>CON2300</b> | C575-20MA  | 20           | 500  | £146.00 |
| <b>CON2302</b> | C575-40MA  | 40           | 500  | £180.00 |
| <b>CON2304</b> | C575-60MA  | 60           | 500  | £170.00 |
| <b>CON2306</b> | C575-90MA  | 90           | 400  | £171.00 |
| <b>CON2308</b> | C575-120MA | 120          | 300  | £148.00 |

## Simport® Scientific SecuriTainer™ III Containers



- Chemically resistant and shatterproof
- Supplied with a leakproof, polyethylene screw cap, particularly important when transporting hazardous material
- Ideal for transport and storage of urine, sputum and most liquids or particulate samples
- Stackable
- Transparent
- Shatter resistant which eliminates issues of leakage and evaporation
- Manufactured from virgin, high clarity polypropylene
- Functional temperature range: -90 to 100°C
- Non sterile



| Code           | Alt Ref   | Capacity, mL | Pack | Price   |
|----------------|-----------|--------------|------|---------|
| <b>CON2310</b> | C577-20W  | 20           | 500  | £131.00 |
| <b>CON2322</b> | C577-40W  | 40           | 500  | £168.00 |
| <b>CON2324</b> | C577-60W  | 60           | 500  | £168.00 |
| <b>CON2326</b> | C577-90W  | 90           | 400  | £170.00 |
| <b>CON2319</b> | C577-120W | 120          | 300  | £184.00 |

## Azlon® Screw Cap Containers



- Wide neck with linerless screw caps
- Polypropylene
- Good chemical resistance
- Ideal for storage of samples, powders, etc
- Autoclavable

| Code           | Alt Ref  | Capacity, mL | Dims, Ø x h, mm | Pack | Price   |
|----------------|----------|--------------|-----------------|------|---------|
| <b>JAR2004</b> | JWP0030P | 30           | 43 x 43         | 72   | £47.10  |
| <b>JAR2006</b> | JWP0060P | 60           | 53 x 47         | 48   | £40.20  |
| <b>JAR2010</b> | JWP0125P | 125          | 70 x 60         | 36   | £40.30  |
| <b>JAR2016</b> | JWP0250P | 250          | 89 x 66         | 36   | £81.00  |
| <b>JAR2022</b> | JWP0500P | 500          | 89 x 97         | 24   | £58.00  |
| <b>JAR2026</b> | JWP0950P | 950          | 120 x 100       | 24   | £151.00 |
| <b>JAR2028</b> | JWP1200P | 1200         | 120 x 132       | 24   | £140.00 |

