

Simport® Scientific Cryovials® Continued



Internal Thread and Silicone Washer Seal


Silicon seal fitted inside the cap ensures a positive seal at any temperature.



Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Shape	Pack	Price
CRY2208	T311-1	1.2	41 x 12.5	Self standing	1000	£320.00
CRY2210	T311-2	2	49 x 12.5	Self standing	1000	£285.00
CRY2212	T311-3	2	48 x 12.5	Round bottom	1000	£322.00
CRY2214	T311-4	4	70 x 12.5	Round bottom	1000	£431.00
CRY2216	T311-4A	4	72 x 12.5	Self standing	1000	£378.00
CRY2218	T311-5	5	90 x 12.5	Round bottom	1000	£391.00

External Thread and Silicone Washer Seal


Silicon seal fitted inside the cap ensures a positive seal at any temperature.



Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Shape	Pack	Price
CRY2220	T310-1A	1.2	42 x 12.5	Self standing	1000	£391.00
CRY2222	T310-2	2	47 x 12.5	Round bottom	1000	£403.00
CRY2224	T310-2A	2	49 x 12.5	Self standing	1000	£348.00
CRY2226	T310-3A	3	71 x 12.5	Self standing	1000	£399.00
CRY2228	T310-4A	4	77 x 12.5	Self standing	1000	£435.00
CRY2230	T310-5A	5	91 x 12.5	Self standing	1000	£504.00
CRY2232	T310-10A	10	84 x 17	Self standing	500	£260.00

External Thread and Lip Seal Design

A specially designed lip inside the cap ensures a leak proof seal even at very low temperatures.



Code	Alt Ref	Capacity, mL	Dims, h x Ø, mm	Shape	Pack	Price
CRY2246	T309-1A	1.2	42 x 12.5	Self standing	1000	£288.00
CRY2248	T309-2	2	47 x 12.5	Round bottom	1000	£233.00
CRY2250	T309-2A	2	49 x 12.5	Self standing	1000	£262.00
CRY2252	T309-3A	3	71 x 12.5	Self standing	1000	£264.00
CRY2254	T309-4A	4	77 x 12.5	Self standing	1000	£290.00
CRY2256	T309-5A	5	91 x 12.5	Self standing	1000	£319.00