

## SLS Select Autoclave Bags, HDPE



- High density polyethylene bags for general disposal
- Autoclavable up to 121°C
- Sizes quoted are approximate

Code	Dims, w x l, mm	Pack	Price
AUT1400	670 x 305	200	£39.40
AUT1402	625 x 406	200	£49.50
AUT1404	840 x 635	100	£48.30



A

## SLS Select Autoclave Bags, PP



- Transparent polypropylene disposal bags
- Autoclavable to 160°C
- Sizes quoted are approximate

Code	Dims, w x l, mm	Pack	Price
AUT1410	634 x 413	500	£138.00
AUT1415	766 x 608	300	£246.00



A

## SLS Select Autoclave Bags, Paper



- Designed for safe, effective packing of instruments, dressings, and pipettes prior to steam sterilisation
- Manufactured from medical grade paper, adhesive and inks
- Designed with a stepped top for easy filling
- Available with a heat seal band and steam indicator or plain top without indicators
- Serrated top to prevent cuts
- Printed with Class 1 process indicators compliant with ISO 11140-1:2014
- Manufactured under ISO 9001:2015 and ISO 13485:2016 quality system
- Conforms to ISO 11607-1:2019 and BS EN 868-4:2017



A

### Heat Seal Closure

Code	Dims, w x d x l, mm	Pack	Price
AUT1432	90 x 50 x 170	1000	£65.00
AUT1434	90 x 50 x 250	1000	£97.00
AUT1435	110 x 30 x 190	1000	£78.00
AUT1440	125 x 30 x 170	1000	£87.00
AUT1442	125 x 50 x 250	1000	£128.00
AUT1444	140 x 75 x 250	1000	£129.00
AUT1445	140 x 50 x 330	500	£77.00
AUT1448	180 x 95 x 380	500	£112.00
AUT1450	190 x 65 x 330	500	£97.00
AUT1451	250 x 100 x 380	250	£75.00
AUT1454	380 x 125 x 610	250	£150.00

### Plain Seal Closure

Code	Dims, w x d x l, mm	Pack	Price
AUT1502	90 x 50 x 170	1000	£71.00
AUT1504	90 x 50 x 250	2000	£95.00
AUT1508	110 x 30 x 190	1000	£76.00
AUT1512	125 x 30 x 170	1000	£81.00
AUT1514	125 x 50 x 250	1000	£123.00
AUT1516	140 x 75 x 250	1000	£124.00
AUT1517	140 x 50 x 330	500	£76.00
AUT1520	180 x 95 x 380	500	£109.00
AUT1521	180 x 75 x 530	250	£79.00
AUT1523	250 x 100 x 380	250	£72.00
AUT1526	380 x 125 x 610	250	£160.00