

CLEARLine® CryoGen® Tubes/2D Barcode Continued

ClearLine



External Thread

Supplied in packs of 500 with inner packs of 50

Code	Alt Ref	Volume, mL	Price
CRY1182	CL2ARBEPS2D	2	£498.00
CRY1184	CL5ARBEPS2D	5	£498.00




Corning® Internal Thread Cryogenic Vials

Corning

- Manufactured from polypropylene to withstand temperatures down to -196°C
- Black graduations with large marking spot
- Vials have a silicone washer for a secure seal
- Vials may be colour coded with inserts
- Free foam rack with each case
- Orange cap
- RNase/DNase free
- Sterile
- Non pyrogenic

Code	Capacity, mL	Bottom	Pack	Price
430487	1.2	Conical	500	£310.00
430488	2	Self standing, round	500	£310.00
431386	2	Self standing, round	250	£159.00
430489	2	Round	500	£310.00
430492	5	Round	500	£395.00



Thermo Scientific™ Mr. Frosty™ Container

Thermo
SCIENTIFIC

- Freeze cells in tubes from 1 to 2mL
- Designed to achieve a rate of cooling very close to -1°C/minute, the optimal rate for cell preservation
- Height saves precious freezer space during the cooling process
- Screw top lid secures samples inside in the event of accidental container tipping
- Eliminates unreliable, improvised foam container devices
- Easy to use, requires only 100% isopropyl alcohol and a mechanical freezer
- No contact with the alcohol: no wicking, no label removal and no compromise of tube hand printing

Code	Alt Ref	Capacity	Dims, h x Ø, mm	Price
CRY8400	5100-0001	18 x 1 to 2mL tubes	86 x 117	£135.00