

Thermo Scientific™ Benchtop Coolers

Thermo
SCIENTIFIC

- Maintain freezer temperature during frost free cycles, power failures and accidental door openings
- Protect enzymes, cells, reagents and solutions on the laboratory bench
- Printed grid on one-way keyed lid for proper orientation and sample identification
- Moulded of impact-resistant polycarbonate with locking handle for long term use
- Suitable for use to -20°C
- Stackable to save space




Code	Alt Ref	Description	Capacity	Array	Dims, w x d x h, mm	Holds <15°C, Hours	Price
CRY8374	5115-0012	Non filled clear lid	12 x 0.2 - 0.5mL tubes	3 x 4	108 x 151 x 125	1	£195.00
CRY8376	5115-0032	Gel filled white lid	32 x 0.2 - 0.5mL tubes	4 x 8	157 x 243 x 146	2	£377.00

Thermo Scientific™ Cryoboxes

Thermo
SCIENTIFIC

- Manufactured from polycarbonate
- Durable, all plastic cryoboxes for short or long term storage
- Keyed lids on plastic cryoboxes to prevent misalignment
- Economical, yet durable, alternative to card or expensive stainless steel
- Autoclavable for decontamination
- Usable at temperature range from -196 to 121°C
- Printed grid on plastic lids for easy inventory management




Code	Alt Ref	Capacity	Array	Dims, w x d x h, mm	Pack	Price
CRY8300	5025-0505	25 x 1.2 & 2mL vials	5 x 5	76 x 76 x 51	8	£78.00
CRY8302	5026-0909	81 x 1.2 & 2mL vials	9 x 9	133 x 133 x 52	4	£63.00
CRY8304	5027-0909	81 x 5mL vials	9 x 9	133 x 133 x 95	4	£78.00
CRY8305	5026-1010	100 x 1 - 1.5mL vials	10 x 10	133 x 133 x 52	10	£162.00

Thermo Scientific™ Nalgene™ Microcentrifuge Tube Storage Boxes

Thermo
SCIENTIFIC

- Bottom constructed of white polycarbonate
- Transparent polycarbonate lid with printed grid for inventory purposes
- High-temperature polyurethane foam insert securely holds microcentrifuge tubes
- Suitable for use from -135 to 121°C
- Lid accepts marker for box and tube identification and is keyed to orient front and back of box
- Dims, w x d x h: 133 x 133 x 51mm

Supplied in packs of 4



Code	Alt Ref	Capacity	Array	Price
CRY8308	5055-5005	81 x 0.5mL tubes	9 x 9	£83.00
CRY8310	5055-5015	64 x 1.5mL tubes	8 x 8	£83.00