




Thermo Scientific™ Undivided Cryobox

Thermo
SCIENTIFIC

- Polycarbonate box
- Suitable for standard sized cryobox racks
- Suitable for mechanical freezer or vapour phase liquid nitrogen storage

Code	Alt Ref	Array	Pack	Price
CRY8306	5050-0001	Undivided	4	£63.00



Thermo Scientific™ Coloured Cryoboxes

Thermo
SCIENTIFIC

- Durable, polycarbonate cryoboxes for short or long term storage
- Keyed lids prevent misalignment
- Standard 9 x 9 configuration with coloured grid for ease of identification
- Compatible with 1 - 2mL vials
- Suitable for use between -196 to 121°C
- Printed grid on plastic lids for easy inventory management
- Autoclavable for decontamination
- Dims, L x w x h: 133 x 133 x 51mm

Code	Alt Ref	Grid Colour	Pack	Price
CRY8280	867013-0240	Red	4	£69.00
CRY8282	867013-0241	Yellow	4	£65.00
CRY8284	867013-0242	Green	4	£65.00
CRY8286	867013-0243	Blue	4	£65.00
CRY8288	867013-0244	Grey	4	£65.00
CRY8289	867013-0245	Assorted (white, red, yellow, green, blue, grey)	24	£389.00



Thermo Scientific™ Nunc™ CryoTube™ Rack

Thermo
SCIENTIFIC

- Autoclavable
- Allows one hand operation when used with starfoot tubes
- Material: polyphenylene oxide

Code	Capacity	Dims, w x d, mm	Price
376589K	4 x 10 tubes	202 x 102	£35.00

Liquid handling... Down to the last drop!

From pipettes to dispensers, tips to pipetting stations,
SLS is your No.1 source for liquid handling.

See pages 745–815



**SCIENTIFIC
LABORATORY
SUPPLIES**

