


Bravo Storage/Cryo Boxes Continued

Dividers

- Height: 30mm


Code	Array	Tube Ø, mm	Price
FRE2220	14 x 14	7	£5.60
FRE2222	13 x 13	8	£5.10
FRE2224	12 x 12	9	£7.10
FRE2226	10 x 10	12	£2.90
FRE2228	9 x 9	13	£2.30
FRE2230	8 x 8	14	£2.20
FRE2232	7 x 7	16	£2.30
FRE2234	6 x 6	19	£2.80
FRE2238	4 x 4	30	£2.90
FRE2236	5 x 5	23	£2.30



CLEARLine® Polycarbonate Cryoboxes

ClearLine

- For storing at temperatures as low as -196°C
- Autoclavable at 121°C for 20 minutes
- Numerical grid printed on the lid
- Frosted area on both lid and bottom for an easy writing
- Removal of vials facilitated by a vial picker supplied with each pack of storage boxes
- Polycarbonate



Code	Alt Ref	Capacity	Dims, w x d x h, mm	Colour	Pack	Price
CRY1114	BSM58025B/AZ	25 x 1 - 2mL tubes	76 x 76 x 52	White with blue grid	8	£38.60
CRY1116	BSM58025B/RO	25 x 1 - 2mL tubes	76 x 76 x 52	White with red grid	8	£35.40
CRY1118	BSM58025B/V	25 x 1 - 2mL tubes	76 x 76 x 52	White with green grid	8	£38.90
CRY1120	BSM58025B/G	25 x 1 - 2mL tubes	76 x 76 x 52	White with yellow grid	8	£35.40
CRY1122	BSM58061B/AZ	81 x 1 - 2mL tubes	132 x 132 x 52	White with blue grid	5	£32.00
CRY1124	BSM58061B/RO	81 x 1 - 2mL tubes	132 x 132 x 52	White with red grid	5	£29.30
CRY1126	BSM58061B/V	81 x 1 - 2mL tubes	132 x 132 x 52	White with green grid	5	£29.30
CRY1128	BSM58061B/G	81 x 1 - 2mL tubes	132 x 132 x 52	White with yellow grid	5	£29.30
CRY1130	BSM58071B/AZ	81 x 3 - 5mL tubes	132 x 132 x 94	White with blue grid	4	£32.10
CRY1132	BSM58071B/RO	81 x 3 - 5mL tubes	132 x 132 x 94	White with red grid	4	£29.40
CRY1134	BSM58071B/V	81 x 3 - 5mL tubes	132 x 132 x 94	White with green grid	4	£32.10
CRY1136	BSM58071B/G	81 x 3 - 5mL tubes	132 x 132 x 94	White with yellow grid	4	£29.40
CRY1138	BSM58100/B	100 x 1 - 2mL tubes	132 x 132 x 52	White	5	£40.90
CRY1140	BSM58100/AZ	100 x 1 - 2mL tubes	132 x 132 x 52	Blue	5	£40.10
CRY1142	BSM58100/RO	100 x 1 - 2mL tubes	132 x 132 x 52	Red	5	£36.80
CRY1144	BSM58100/V	100 x 1 - 2mL tubes	132 x 132 x 52	Green	5	£36.80
CRY1146	BSM58100/G	100 x 1 - 2mL tubes	132 x 132 x 52	Yellow	5	£36.80