




CLEARLine® 1D Polycarbonate Cryoboxes

ClearLine

- Linear barcode 'code 128' lasered on the rack side
- Suitable for automatic storage system
- Numerical grid serigraphed both on lid and bottom for an easy positioning
- Frosted area on both lid and bottom for an easy writing
- Coloured grid
- For storing at temperatures as low as -196°C
- Autoclavable at 121°C for 20 minutes
- Removal of vials facilitated by a vial picker supplied with each pack of storage boxes
- Polycarbonate

Code	Alt Ref	Capacity	Dims, w x d x h, mm	Grid Colour	Pack	Price
CRY1148	BSM58025B1D/AZ	25 x 1 - 2mL tubes	76 x 76 x 52	Light blue	8	£82.30
CRY1150	BSM58025B1D/RO	25 x 1 - 2mL tubes	76 x 76 x 52	Red	8	£76.20
CRY1152	BSM58025B1D/V	25 x 1 - 2mL tubes	76 x 76 x 52	Green	8	£82.30
CRY1154	BSM58025B1D/G	25 x 1 - 2mL tubes	76 x 76 x 52	Yellow	8	£82.30
CRY1156	BSM58061B1D/AZ	81 x 1 - 2mL tubes	132 x 132 x 52	Light blue	5	£73.30
CRY1158	BSM58061B1D/RO	81 x 1 - 2mL tubes	132 x 132 x 52	Red	5	£73.30
CRY1160	BSM58061B1D/V	81 x 1 - 2mL tubes	132 x 132 x 52	Green	5	£80.00
CRY1162	BSM58061B1D/G	81 x 1 - 2mL tubes	132 x 132 x 52	Yellow	5	£73.30
CRY1164	BSM58071B1D/AZ	81 x 3 - 5mL tubes	132 x 132 x 94	Light blue	4	£64.00
CRY1166	BSM58071B1D/RO	81 x 3 - 5mL tubes	132 x 132 x 94	Red	4	£58.60
CRY1168	BSM58071B1D/V	81 x 3 - 5mL tubes	132 x 132 x 94	Green	4	£58.60
CRY1170	BSM58071B1D/G	81 x 3 - 5mL tubes	132 x 132 x 94	Yellow	4	£58.60



CLEARLine® 2D Polycarbonate Cryoboxes

ClearLine

- Automatic identification of the box and of its positioning through 2D datamatrix code, on the box bottom
- Barcode "code 128" lasered on rack side
- Numerical grid on lid
- Frosted writing area, white bottom and coloured grid
- For storing at temperatures as low as -196°C
- Autoclavable at 121°C for 20 minutes
- Vial picker supplied with each inner pack
- Polycarbonate

Code	Alt Ref	Capacity	Dims, w x d x h, mm	Grid Colour	Pack	Price
CRY1198	BSM580612D/G	81 x 2mL tubes	132 x 132 x 52	Yellow	5	£80.00
CRY1192	BSM580612D/AZ	81 x 2mL tubes	132 x 132 x 52	Light blue	5	£73.30
CRY1194	BSM580612D/RO	81 x 2mL tubes	132 x 132 x 52	Red	5	£73.30
CRY1196	BSM580612D/V	81 x 2mL tubes	132 x 132 x 52	Green	5	£73.30
CRY1200	BSM580712D/AZ	81 x 5mL tubes	132 x 132 x 94	Light blue	4	£58.60
CRY1202	BSM580712D/RO	81 x 5mL tubes	132 x 132 x 94	Red	4	£58.60
CRY1204	BSM580712D/V	81 x 5mL tubes	132 x 132 x 94	Green	4	£58.60
CRY1206	BSM580712D/G	81 x 5mL tubes	132 x 132 x 94	Yellow	4	£58.60
CRY1208	BSM581002D/B	100 x 2mL tubes	132 x 132 x 52	White	5	£73.30
CRY1210	BSM581052D/B	200 x 5mL tubes	132 x 132 x 94	White	4	£58.60