

SLS Lab Pro Statebourne Biorack Storage Refrigerators Continued



Code	MOL3146	MOL3182	MOL3183	MOL3184
LN ₂ Liquid Phase Capacity, L	35	85	131	160
LN ₂ Vapour Phase Capacity, L	14	25	25	25
Evaporation Rate, L/day	0.3	0.63	0.63	0.63
Static Holding Time, days	117	135	208	254
Working Time*, days	73	84	130	158
Weight (empty), kg	19	42	46	64
Weight (full), kg	47	110	152	193
Overall Dims, h x Ø, mm	710 x 480	738 x 646	913 x 646	1027 x 648
Neck Ø, mm	120	216	216	216
Vial Capacity, 2mL	750	3000	4800	6000

* Actual working time will depend on applications

Code	Alt Ref	Model	Supplied with	Price
MOL3146	9902132	Biorack 750	5 x racks, 30 x 5 x 5 cryoboxes	£2933.00
MOL3182	9902149	Biorack 3000	6 x racks, 30 x 10 x 10 cryoboxes	£4425.00
MOL3183	9902154	Biorack 4800	6 x racks, 48 x 10 x 10 cryoboxes	£5545.00
MOL3184	9902144	Biorack 6000	6 x racks, 60 x 10 x 10 cryoboxes	£5992.00

Accessories

Code	Alt Ref	Description	For use with	Price
CRY3114	9902112	Roller base, 110mm (h)	Biorack 750	£358.00
CRY3116	9902147	Roller base, 110mm (h)	Biorack 3000, 4800, 6000	£528.00

SLS Lab Pro Statebourne Biostor Wide Neck Refrigerators



A range of wide neck refrigerators for laboratory process, bench top and industrial freezing applications.

- Liquid or vapour phase storage
- Low cost
- High thermal efficiency
- Supplied with built in liquid nitrogen fill and level sensor tubes to allow the addition of an autofill kit using a suitable Cryostor as a reservoir
- Level sensors can be set to control liquid levels giving either liquid or vapour phase storage
- Temperatures as low as -185°C can be achieved at warmest point of sample storage



5
YEARS

5 year warranty on vacuum

Code	MOL3120	MOL3103	CRY3105
LN ₂ Capacity, L	66	220	380
Weight (empty), kg	28	125	162
Weight (full), kg	81	303	469
Internal Dims, Ø x h, mm	406 x 513	617 x 745	788 x 835
Overall Dims, Ø x h, mm	460 x 777	750 x 1141	900 x 1193
Vial Capacity, 2mL	5292	12300	24050
Straw Capacity, 0.25cc	56580	170560	401800
Straw Capacity, 0.5cc	25185	75920	178850

Code	Alt Ref	Model	Price
MOL3120	9915120	Biostor 5	£2790.00
CRY3103	9916027	Biostor 12	£6274.00
CRY3105	9916030	Biostor 24	£10292.00