

SLS Lab Pro Statebourne LN<sub>2</sub> Storage Bio Series Continued

## Long Holding Time Refrigerators Continued

Code	Alt Ref	Model	Price
CRY3000	99-02015	Bio 2	£775.00
CRY3002	99-02030	Bio 3	£857.00
CRY3004	99-02050	Bio 8	£1010.00
CRY3006	99-02060	Bio 10	£1050.00
CRY3010	99-02080	Bio 20	£1085.00
CRY3016	99-02120	Bio 35	£1456.00

## Accessories

Code	Alt Ref	Description	For use with	Pack	Price
CRY3017	97-1065	Tote bag	Bio 2	Each	£78.00
CRY3019	97-1011	Tote bag	Bio 3	Each	£129.00
MOL3124	9702000	Cryoguard alarm	Bio 8, 10, 20, 34, 35, 36	Each	£860.00
CRY1104	99-02082	Roller base, 110mm (h)	Bio 8, 12, 20	Each	£362.00
CRY3114	9902112	Roller base, 110mm (h)	Bio 21, 22, 34, 35, 36	Each	£358.00
CRY3021	5801002	Ampoule canes	Bio 20, 34, 35, 36	100	£264.00

SLS Lab Pro Statebourne Cryoguard LN<sub>2</sub> Low Level Alarm

For use with liquid nitrogen refrigerators Bio, Biorack and Biostor. Supplied with control box and adjustable probe fitted in a protective tube specially mounted in the neck plug of the vessel. When the level of liquid nitrogen falls below the set position of the probe, both visual and audible alarms are activated.

Code	Alt Ref	Power Supply	Price
MOL9012	8602622	Mains or battery powered	£860.00



## Accessories

Code	Alt Ref	Description	Price
MOL9020	9702061	Alarm plug BR24K/30K/48K/54K/60K	£307.00
MOL9022	9702155	Alarm plug BR750	£189.00

## KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty