

SLS Select Sterilisation Pouches



- Self seal pouches with a self adhesive strip which creates a seal suitable for withstanding autoclave or ethylene oxide gas sterilisation
- Colour change indicator arrows

Supplied in packs of 200

Code	Dims, w x l, mm	Price
AUT1550	133 x 90	£39.00
AUT1551	203 x 90	£34.00
AUT1552	230 x 90	£35.00
AUT1554	254 x 133	£57.00
AUT1556	330 x 190	£65.00
AUT1558	381 x 305	£115.00



A

SLS Select Sterilisation Bag Reels



- Pouches for autoclave or ethylene oxide gas sterilisation
- Meets the requirements of BS 6871
- Colour change indicator arrows
- Supplied as a reel to be cut as required

Code	Dims, mm x m	Price
AUT1562	75 x 100	£100.00
AUT1564	100 x 100	£107.00
AUT1566	150 x 100	£129.00
AUT1568	200 x 100	£128.00
AUT1570	250 x 100	£166.00
AUT1572	300 x 100	£186.00



A

Heat Sealers

- Designed primarily for sealing polythene and polypropylene film
- Sealing is only initiated when the jaws are closed on the film, with no warm up time required
- Variable electronic timer allows the required seal time to be set relevant to the thickness of the film, ensuring a hermetic seal every time
- LED indicates the duration of the cycle
- Cutter version is available where there is a need to trim bags, the integral blade trimmer is manually operated and cuts the film 10mm behind the seal line
- Simple to operate, portable and compact
- Supplied with three pin plug, ready for use
- Seal width: 2.5mm



Code	Model	Seal Length, mm	Cut Length, mm	Power, W	Dims, w x d x h, mm	Price
BAG1050	205	200	—	280	150 x 80 x 320	£181.00
BAG1052	205CUTTER	200	200	280	180 x 80 x 370	£245.00
BAG1054	305	300	—	350	180 x 80 x 440	£263.00
BAG1056	305CUTTER	300	300	350	190 x 90 x 470	£364.00
BAG1058	405	400	—	450	170 x 90 x 540	£304.00

Polythene Layflat Tubing (LFT)

- Use with above heat sealer
- Cut to size

Code	Dims, mm x m	Price
BAG1072	150 x 150	£61.00
BAG1074	250 x 150	£115.00
BAG1076	305 x 150	£125.00