



Epredia™ CytoSpin™ 4 Cytocentrifuge



Gain all the advantages of the ultimate thin layer cell preparation system with this reliable benchtop centrifuge. It provides economical thin layer preparations from any liquid matrix, especially hypocellular fluids such as spinal fluid and urine. Processes up to twelve specimens at one time.

- Accepts all protocols from CytoSpin™ 1, 2 and 3
- Allows for one handed opening and closing, only operates when the lid is locked into position and lid remains locked during rotation
- Sealed head can be viewed through the polycarbonate window during operation
- Designed for easy disinfection
- Meets IEC 61010 safety standards
- Sealed head can be loaded inside a microbiological cabinet for safety
- Audible alarm sounds when it is out of balance or outside the speed tolerance
- Specimen safety alarms sound at one minute intervals to remind users to remove specimens to prevent air drying
- Cells deposited onto a clearly defined area of a glass slide, residual fluid is absorbed into the chamber's filter card
- In the Load or Stop position, the tilting feature reduces cell loss by preventing residual fluid from coming into contact with the prepared slide
- Pull out program card logs up to nine programs for fast, convenient retrieval
- Protected program memory stores up to 23 routines for instant recall
- Easy touch program keypad allows for intuitive programming
- Window displays programmed spinning speeds from 200 to 2000RPM and programmed/remaining time
- Three acceleration control buttons allow the operator the flexibility to protect fragile samples

Code	Alt Ref	Dims, w x d x h*, mm	Weight, kg	Price
HIS4030	A78300002	405 x 620 x 240	12	£9556.00

* With closed lid, height of 625mm with open lid

Accessories

Starter kit includes 50 single white cytofunnels, 25 double cytofunnels, cytoslides, and 12 cytoclips.

Code	Alt Ref	Description	Price
HIS4032	3120110	Cytospin 4 clipped funnel starter kit	£1083.00



Simport® Scientific CytoSep™ Cytology Funnels



- Fully compatible with the Shandon CytoSpin® centrifuge
- Can be used to deposit a thin layer of cells in a clearly defined area of a microscope slide
- Filter card absorbs any excess fluid
- Pre-attached filter cards for consistent, reliable results - no alignment necessary
- Packaged with closure caps to seal in the specimen for added protection
- For sample volumes of up to 0.5mL
- Provides a cell deposition area of 6mm (28mm²)
- Lower risk of contamination than reusable funnels
- Double funnels allow for two samples to be run simultaneously on one slide
- Can be used with all stainless steel slide clips

Code	Alt Ref	Description	Pack	Price
CYT4000	M964-10FW	Single cytology funnel with white filter card and cap	50	£70.00
CYT4002	M964-10FW	Single cytology funnel with white filter card and cap	500	£567.00
CYT4005	M964-20FW	Double cytology funnel with white filter card and cap	50	£96.00
CYT4006	M964-20FW	Double cytology funnel with white filter card and cap	500	£788.00