

BRAND Dispensette® S Bottle Top Dispensers Continued



Dispensette® S Organic Fixed-Volume

Code	Alt Ref	Recirculation Valve	Volume, mL	Price
DIS1756	4630 230	No	5	£332.00
DIS1758	4630 240	No	10	£332.00
DIS1760	4630 231	Yes	5	£382.00
DIS1762	4630 241	Yes	10	£382.00



Dispensette® S Trace Analysis Analogue-Adjustable

- Provides outstanding performance for precise volume dispensing of high-purity media for trace analysis
- Plastics in contact with media consist of high-purity materials such as PTFE, ETFE, FEP and PFA
- Purest sapphire is used for replaceable valves



Code	Alt Ref	Valve Spring	Recirculation Valve	Volume, mL	Graduations, mL	Price
DIS2008	4640 040	Platinum-iridium	No	1 - 10	0.2	£715.00
DIS2010	4640 041	Platinum-iridium	Yes	1 - 10	0.2	£760.00
DIS2012	4640 240	Tantalum	No	1 - 10	0.2	£715.00
DIS2014	4640 241	Tantalum	Yes	1 - 10	0.2	£760.00

Accessories

Code	Alt Ref	Description	Price
DIS2016	708022	Discharge tubes pt-lr w/o re-circulation valve	£15.70
DIS2018	708122	Discharge tubes pt-lr w re-circulation valve	£94.00
DIS2020	708024	Discharge tubes ta valve spring w/o re-circulation valve	£15.70
DIS2022	708124	Discharge tubes ta valve spring w re-circulation valve	£94.00
DIS2024	708035	Dispensing cartridge with safety ring	£423.00
DIS2026	708210	Telescoping filling tubes 70 - 140mm	£12.50
DIS2028	708212	Telescoping filling tubes 125 - 240mm	£17.70
DIS2030	708214	Telescoping filling tubes 195 - 350mm	£25.00
DIS2032	708216	Telescoping filling tubes 250 - 480mm	£19.80