


INTEGRA Biosciences VACUBOY Hand Operator

INTEGRA

A convenient tool that converts a vacuum source into a user friendly system for aspirating liquids from virtually any laboratory container - at the touch of a button and in a controlled way. The hand piece, the core element of the VACUBOY, is simply connected to any available vacuum source using silicone tubing and fitted with one of the different aspiration adapters. The perfect match is obtained when the VACUBOY is connected to a VACUSAFE aspiration system.

- Large choice of adapters allows VACUBOY to handle liquids from all tubes, dishes, bottles, T-flasks or multi-well plates, suiting their form and size
- Adapters for aspiration with Pasteur pipettes, serological pipettes, disposable pipette tips or stainless steel needles
- Multi-channel adapters for fast emptying of multi-well plates
- Vacuum regulation by a touch-sensitive button for controlled and smooth aspiration of liquids
- Ergonomic, perfectly suited for long working, productive or repetitive aspiration tasks
- Anti-drip system prevents liquid waste leakage from the tip or pipette
- Autoclavable and sterilisable for maximum safety in the lab



Code	Description	Price
I155510	VACUBOY hand operator	£158.00
I155500	VACUBOY hand operator, stand for hand operator, rubber adapter for Pasteur pipettes, 1-channel stainless steel adapter 40mm, 8-channel stainless steel adapter 40mm, 1-channel adapter for disposable tips/GripTips (with locking), 1-channel adapter with ejector for disposable tips (with locking), 8-channel adapter with ejector for disposable tips (with locking), silicone tubing (180cm), silicone grease	£267.00

Accessories and Spare Parts

Code	Description	Price
I155501	Stand for VACUBOY	£83.00
I155502	1-channel stainless steel adapter, 40mm	£13.20
I155503	8-channel stainless steel adapter, 40mm	£70.00
I155504	1-channel adapter for disposable tips/GRIPTIPS (with locking)	£16.80
I155505	Rubber adapter for Pasteur pipettes	£9.80
I155520	8-channel adapter with ejector for disposable tips (with locking)	£120.00
I155522	1-channel stainless steel adapter, 150mm	£26.70
I155524	4-channel stainless steel adapter, 40mm	£59.00
I155525	1-channel stainless steel adapter, 280mm	£31.70
I155526	1-channel adapter with ejector for disposable tips (with locking)	£28.90

KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty