



Whirl-Pak® Sampling Bags

- Blended, virgin polyethylene film for exceptional clarity and superior strength
- Sterilised after manufacturing for increased sample integrity
- Leakproof closure: whirling or folding the bag over three times and securing with the closure tabs creates a leakproof sampling container
- Continuous perimeter seal for maximum protection against leaks
- Puncture proof closure tabs for easy access and no leak closure with tape that extends beyond wire to eliminate sharp ends and provide a sturdy grip, wire tabs
- Non slip, textured Tyvek opening pull tabs for easy opening
- Can be used at temperatures up to 82°C, and can be frozen to any temperature, even liquid nitrogen temperatures down to -210°C
- Capacity is measured when the bag is closed and tab is folded over three times
- Volume and dimensions are approximate

Standard Bags

Code	Alt Ref	Description	Capacity, mL	Dims, w x l, mm	Pack	Price
BAG2502	B01009	Round wire	58	75 x 125	500	£82.00
BAG2506	B00679	Round wire	118	75 x 185	500	£92.00
BAG2510	B00992	Round wire	207	95 x 180	500	£115.00
BAG2512	B01018	Round wire	384	130 x 190	500	£126.00
BAG2514	B00736	Round wire	532	115 x 230	500	£116.00
BAG2518	B01020	Round wire	710	150 x 230	500	£166.00
BAG2522	B01063	Round wire	710	150 x 230	500	£205.00
BAG2528	B00990	Round wire	798	125 x 305	500	£155.00
BAG2532	B00994	Round wire	1065	125 x 380	500	£239.00
BAG2534	B01027	Round wire	1242	150 x 380	500	£214.00
BAG2536	B01323	Round wire	2041	190 x 380	500	£287.00



Write On Bags

- Write on strip makes labeling quick and easy for clear sample identification, organisation and traceability
- Flat wire provides a thicker tab, resulting in a heavier closure

Code	Alt Ref	Description	Capacity, mL	Dims, w x l, mm	Pack	Price
BAG2500	B01067	Round wire	29	65 x 125	500	£100.00
BAG2504	B01064	Round wire	58	75 x 125	500	£92.00
BAG2508	B01062	Round wire	118	75 x 185	500	£103.00
BAG2516	B01065	Round wire	532	115 x 230	500	£129.00
BAG2520	B01297	Round wire	710	150 x 230	500	£172.00
BAG2524	B01196	Round wire	710	150 x 230	500	£214.00
BAG2530	B01195	Round wire	1627	190 x 300	500	£298.00
BAG2538	B01445	1 flat, 1 round wire	2721	254 x 380	250	£309.00
BAG2540	B01446	1 flat, 1 round wire	3637	254 x 508	250	£309.00
BAG2542	B01447	1 flat, 1 round wire	5441	380 x 508	100	£192.00
BAG2550	B01339	Flat wire	118	70 x 170	500	£145.00
BAG2552	B01341	Flat wire	532	115 x 230	500	£162.00