

Hanna edge® Hybrid Multiparameter Series Continued



Model	pH Range	mV Range, mV	Temperature, °C	EC Range, $\mu\text{S}/\text{cm}$	TDS Range, mg/L	Salinity	DO Range, ppm
PHM4300	-2 to 16, ± 0.01	$\pm 1000, 0.1$	-20 to 120, ± 0.5	—	—	—	—
PHM0342	—	—		0.00 to 29.99, $\pm 1\%$ of reading	0.00 to 14.99, $\pm 1\%$ of reading	0.0 to 400.0‰ NaCl, $\pm 1\%$ of reading	—
PHM0344	—	—		—	—	—	0.00 to 45.00, $\pm 1.5\%$ reading

Ranges are relevant for all meters but require purchase of extra probes to measure where no value is listed

Code	Alt Ref	Description	Supplied with	Price
PHM4300	HI-2020-02	edge® Hybrid Multiparameter pH meter	pH probe, pH calibration, cleaning sachets	£532.00
PHM0342	HI-2030-02	edge® Hybrid Multiparameter EC/TDS/salinity meter	EC/TDS/salinity probe, EC calibration sachets	£755.00
PHM0344	HI-2040-02	edge® Hybrid Multiparameter DO meter	DO probe, electrolyte solution, 2 x DO membranes	£755.00

Hanna edge® blu Bluetooth pH Meter



Clearing the clutter, edge® blu sees the end of trailing cables by wirelessly connecting to the pH electrodes with Bluetooth® Smart technology. A hybrid meter that can be used in portable, wall-mount and benchtop configurations. The versatile design of edge® pH enables it to be used as a portable, wall-mount or benchtop meter. edge® pH simplifies measurement, configuration, calibration, diagnostics, logging and transferring data directly to a computer or USB drive.

- Two USB ports, one standard, one micro
- Store up to 1000 data sets including readings, GLP data, date and time
- Exclusive CAL Check™ technology warns you if the electrode bulb is not clean or if the buffers are contaminated during calibration
- Capacitive touch keypad
- Easy to read 14cm LCD with wide viewing angle
- Incredibly thin and lightweight, measures just 12mm thick and weighs just 250g
- pH electrode with Bluetooth® Smart technology



Specifications	
pH Range	Basic mode: -2.00 to 16.00, standard mode: -2.000 to 16.000
pH Resolution	0.001, 0.01
pH Accuracy	$\pm 0.01, \pm 0.002$
Temperature Range, °C	-20.0 to 120.0
Temperature Resolution, °C	0.1
Temperature Accuracy, °C	± 0.5

Code	Alt Ref	Description	Supplied with	Price
PHE1038	HI-2202	edge® blu pH meter	Bench-top stand, electrode holder, bracket for wall-mount, USB cable, power adaptor, instruction manual, HI-11102 pH electrode, pH 7 and 4 buffer solutions, electrode cleaning solution	£568.00

Oxoid and Difco media available!

A wide range of Oxoid and Difco media are available.

Visit the SLS website to view the range.

