



Mettler Toledo InLab® pH Electrodes



The complete InLab® pH sensor portfolio covers the whole scope of applications. These range from routine measurements to specialist applications, be it in chemical, pharmaceutical, life science, food and beverage or a variety of other laboratory or production facilities.

- ARGENTHAL™ reference system ensures the electrolyte will remain completely free of silver ions by including a silver ion trap
- Included wetting cap ensures your pH sensor is kept hydrated and is ready to use at any time
- Different styles of junctions to ensure performance regardless of the sample type
- Easy to use and spill proof cap ensures your sensor will remain clean and secure during storage

Code	Alt Ref	Description	Range, pH	Temp Range, °C	Connector	Material	Junction	Reference Electrolyte	Dims, w x l, mm	Price
PHE0300	51343011	Easy BNC	0 - 14	0 - 80	BNC	PSU	Ceramic	Gel	12 x 120	£125.00
PHE0302	51343010	Easy	0 - 14	0 - 80	S7	PSU	Ceramic	Gel	12 x 120	£97.20
PHE0306	51343101	Expert Pro*	0 - 14	0 - 100	BNC/RCA (cinch)	PEEK	Open	XEROLYT® Polymer	12 x 120	£244.00
PHE0310	51343054	Routine Pro*	0 - 14	0 - 100	MultiPin™	Glass	Ceramic	3M KCl	12 x 120	£239.00
PHE0312	51343056	Routine**	0 - 14	0 - 100	MultiPin™	Glass	Ceramic	3M KCl	12 x 120	£239.00
PHE0314	51343050	Routine	0 - 14	0 - 100	S7	Glass	Ceramic	3M KCl	12 x 120	£161.50
PHE0318	51343105	Expert Pro**	0 - 14	0 - 100	MultiPin™	PEEK	Open	XEROLYT® Polymer	12 x 120	£210.00
PHE0320	51343070	Science	0 - 12	0 - 100	S7	Glass	Movable glass sleeve	3M KCl	12 x 120	£235.00
PHE0322	51343165	Semi-Micro	0 - 12	0 - 100	S7	Glass	Open	XEROLYT® EXTRA Polymer	6 x 100	£244.00
PHE0324	51343160	Micro	0 - 14	0 - 80	S7	Glass	Ceramic	3M KCl	3 x 60	£272.00
PHE0326	51343157	Surface	1 - 11	0 - 50	S7	Glass	Ceramic	3M KCl	12 x 120	£229.00
PHE0328	51343153	Solids	1 - 11	0 - 80	S7	Glass	Open	XEROLYT® EXTRA Polymer	6 x 25	£272.00
PHE0330	51343174	Cool	1 - 11	-30 - 80	S7	Glass	Immovable glass sleeve	FRISCOLYT-B®	12 x 120	£260.00
PHE0332	51343176	Hydrofluoric	1 - 11	0 - 100	S7	Glass	Ceramic	3M KCl	12 x 120	£284.00
PHE0334	51343161	Semi-Micro-L	0 - 14	0 - 100	S7	Glass	Ceramic	3M KCl	6 x 230	£270.00
PHE0338	51343062	Reach Pt1000-425**	0 - 14	0 - 100	MultiPin™	Glass	Ceramic	3M KCl	12 x 425	£306.00
PHE0340	51343061	Reach Pro-425*	0 - 14	0 - 100	MultiPin™	Glass	Ceramic	3M KCl	12 x 400	£306.00
PHE0344	51343195	Mono	0 - 14	0 - 100	S7	Glass	—	—	12 x 120	£119.50
PHM2636	51343164	Flex-Micro	0 - 14	0 - 80	BNC	Epoxy	Porous PTFE	Gel	6 x 180	£360.00

* Temperature probe: NTC 30 kΩ

** Temperature probe: Pt1000

Take a closer look at our microscope range!

Including monocular, binocular and trinocular models.

See pages 623–649



**SCIENTIFIC
LABORATORY
SUPPLIES**

