



Sartorius Secura® Semi-Micro Balances

SARTORIUS

Reduce your risk with integrated and intelligent protection systems. Secura® gives you the security and peace of mind of knowing that you have done everything right. Besides providing highly accurate weighing results and operating convenience, built-in protection systems are included for complete reliability and regulatory compliance – the safe and secure way. Real time guidance prompts for levelling, automatic internal adjustment, and 100% traceable, clear documentation with sample and batch identifiers make your lab work more efficient.

- Integrated protection systems for minimising your risk of obtaining incorrect results
- Real time alarm messaging and interactive user guidance for levelling to prevent inaccurate results due to an unlevelled balance
- SQmin for automatic detection and active monitoring to ensure you are within the acceptable USP operating range for reliable results
- Fully automatic internal adjustment function, isoCAL, for the highest possible accuracy at all times
- Cal Audit Trail for easy, traceable documentation of all calibration and adjustment data
- Safe and clear documentation by sample and batch identifiers for keeping perfect track of weighing data
- GLP/GMP compliant printouts
- Standard password security features for protection against tampering and unauthorised access
- Integrated PC direct function and plug and play connectivity to a Sartorius printer
- Pan dims: Ø 80mm

Code	Alt Ref	Capacity, g	Readability, mg	Price
BAL1158	SECURA125-1S	120	0.01	£6845.00
BAL3874	SECURA225D-1S	120/220	0.01/0.1	£6882.00



OHAUS Pioneer® PX Semi-Micro Balances



Combine essential weighing functionality with excellent performance, at an affordable price and is carefully designed for durability and ease of use.

- High accuracy and repeatability for essential weighing applications
- Cast metal lower housing, sub-pan and stainless steel weighing pan
- Durably constructed for versatile, long term use
- 2nd line display for additional information or guidance, a static removal bar for convenient grounding, and USB connectivity
- Other features include: RS232 interface, integral weigh below hook, removable stainless steel pan, security bracket, software lockout menus, stability indicator, software overload/underload indicators, user selectable environmental settings, user brightness settings, auto dim, auto tare, user selectable operating language (EN, CN, JP, KR)
- Internal calibration
- Pan dims: Ø 80mm

Code	Alt Ref	Capacity, g	Readability, mg	Price
BAL1312	PX125D	52/120	0.01/0.1	£2545.00
BAL1314	PX85	82	0.01	£2640.00
BAL1316	PX225D	82/220	0.01/0.1	£2735.00

Accessories

Code	Alt Ref	Description	Price
BAL1340	80251396	Auxiliary display PAD	£64.00
BAL1342	83034024	Calibrated sinker for liquid density	£305.00
BAL1344	80500525	RS232 cable (9 pin)	£80.70
BAL5146	83021085	USB interface cable	£17.30
BAL5152	80850043	Security device	£157.50
BAL1021	30064202	Printer SF40A	£466.00
BAL5524	80253384	Density kit	£510.00