

2
YEARS

SLS Flowgen Semi Dry BlotTERS



Suitable for use with all types of blotting, semi dry blotters offer rapid transfer times for DNA, RNA and protein blotting typically 15 - 30 minutes.

- Semi dry and therefore economical on buffer consumption
- Corrosion resistant plate electrodes maximise the intensity of the electric field, facilitating Western, Northern and Southern transfer in <30 minutes
- Screw down lid ensures even transfer of up to 4 mini-gels (SD20) of 10mm max thickness
- Uniform heat dispersion
- Rapid transfer times
- Long life electrodes

Code	Description	Capacity, mm	Dims, w x d x h, mm	Price
SD10	Mini	100 x 100	165 x 150 x 55	£476.00
SD20	Midi	200 x 200	260 x 260 x 70	£1020.00

1
YEAR

Cleaver multiSUB Mini One Horizontal Electrophoresis System



The multiSUB Mini One is an all in one horizontal electrophoresis unit. Featuring a built in power supply with voltage options of 35V, 50V and 100V, it is a versatile system suitable for a wide range of applications. A built in timer for runs from 0-99 minutes means you can set up your parameters and leave the system to complete the run automatically without fear of losing bands. The system comes complete with 2 gel casters for wide and mini gels as well as reversible combs for high throughput or high sample volume.

- All in one horizontal electrophoresis system
- In built power supply with 35V, 50V and 100V settings
- Timer function for runs from 0 - 99 minutes
- 2 gel sizes and reversible comb options
- Bromophenol blue migration rate: ~4cm/hr at 50V
- Typical running conditions: 50V, 30 - 60 minutes

Code	Alt Ref	Gel Dims, w x l, mm	Dims, w x d x h, mm	Max. Capacity	Base Buffer Volume, mL	Price
ELE3400	MSMINIONE	110 x 60 & 54 x 60	190 x 130 x 55	22 samples	230	£548.00

Accessories

Code	Alt Ref	Description	Price
ELE3412	MSO-UVS	Gel tray, small, 54 x 60mm	£26.00
ELE3410	MSO-UVL	Gel tray, large, 110 x 60mm	£29.00
ELE3404	MSO-1-5/9DS	Double comb for small gel tray	£29.00
ELE3402	MSO-1-12/22DS	Full length comb for large gel tray	£29.00
ELE3408	MSO-GCS	Gel caster, small	£33.00
ELE3406	MSO-GCL	Gel caster, large	£33.00