



Lonza FlashGel™ Reagents

Lonza

Reagents for use with the Lonza FlashGel™ system.

Code	Alt Ref	Description	Pack	Price
LZ50462	50462	Loading dye	5 x 1mL	£176.00
LZ50475	50475	Quant ladder 100bp (3ng) - 1.5kb (30ng)	250µL	£193.00
LZ50577	50577	RNA marker (50 applications)	50µL	£214.00
LZ57033	57033	DNA marker 50bp - 1.5kb	500µL	£193.00
LZ57034	57034	DNA marker 100bp - 3kb	500µL	£232.00
LZ50473	50473	DNA marker 100bp - 4kb	500µL	£232.00
LZ57050	57050	Recovery kit (9 x FlashGel™ recovery cassettes 1.2% 8+1 well double tier, FlashGel™ recovery buffer, FlashGel™ loading dye, FlashGel™ QuantLadder, visualization glasses)	Each	£908.00
LZ57060	57060	Recovery buffer	2 x 500µL	£163.00
LZ57063	57063	DNA kit (FlashGel™ 1.2% DNA cassettes, loading dye and marker)	Each	£492.00



Lonza FlashGel™ Accessories

Lonza

Accessories for the Lonza FlashGel™ system.

Code	Alt Ref	Description	Price
LZ57040	57040	FlashGel™ camera	£825.00
LZ57068	57068	FlashGel™ power supply	£592.00
LZ57025	57025	FlashGel™ dock	£591.00
LZ57061	57061	FlashGel™ visualisation glasses	£41.10



SLS Flowgen omniDOC Gel Doc Systems

FLOWGEN

The omniDOC offers high performance gel documentation and analysis at a relatively low cost. By providing many of the features used by leading gel documentation brands, but without the added price premium, each omniDOC system presents a simple but sophisticated imaging solution for most laboratories. A high resolution 5MP monochrome camera with slide out UV transilluminator, and optional blue epi-illumination module and white light table, makes the omniDOC suitable for imaging most fluorescent and colorimetric gels, while a USB port requires a cable to connect the dark room assembly to an external PC for control. Imaging applications are made easy by a pre-focused camera that requires little manual adjustment, while simple one click image acquisition and analysis software guides the user through every step of the gel documentation process. A front LED indicator panel reveals at a glance the light source in use, whereas a viewing screen with universal filter and spring-loaded cover facilitates safe and convenient gel inspection.

- Pre-focused 5 megapixel monochrome camera
- Interchangeable 4 position filter wheel with 620nm as standard
- Viewing pane with universal amber filter
- Internal white LED aids gel
- Slide out 312nm transilluminator
- Optional white light table
- Optional blue light modules
- Large 21 x 26cm filter area

Code	White Light Table	Blue Light	Price
ELE2500	No	No	£5907.00
ELE2502	No	Yes	£6573.00
ELE2504	Yes	No	£6240.00
ELE2506	Yes	Yes	£6682.00