


## Sartorius Secura® Analytical Balances

SARTORIUS

Reduce your risk with integrated and intelligent protection systems. Secura® gives you the security and peace of mind of knowing that you have done everything right. Besides providing highly accurate weighing results and operating convenience, built in protection systems are included for complete reliability and regulatory compliance – the safe and secure way. Real time guidance prompts for levelling, automatic internal adjustment, and 100% traceable, clear documentation with sample and batch identifiers make your lab work more efficient.

- Data transfer to a printer or a computer is blocked when most common operating errors occur
- LevelControl, an optoelectronic level monitoring helps prevent incorrect results from an unlevelled balance
- Active monitoring of the USP minimum sample weight requirement
- isoCAL helps prevent incorrect results due to an uncalibrated balance
- CAL Audit Trail provides traceable documentation of all calibration and adjustment data inside the balance
- GLP compliant printouts with time and date stamps with a Sartorius GLP printer YDP30
- Password protection prevents unauthorised changes of weighing settings
- Easy cleaning thanks to innovative housing design
- Optimised draft ring and LevelControl ensure suitability for use in laminar airflow environments
- Pan dims: Ø 90mm



Code	Alt Ref	Capacity, g	Readability, mg	Price
<b>BAL7000</b>	SECURA124-1S	120	0.1	£4189.00
<b>BAL7004</b>	SECURA224-1S	220	0.1	£4513.00

## Sartorius Practum® Analytical Balances

SARTORIUS

Get the best value for your money, without compromises in precision and reliability. Rely on consistent readings and excellent repeatability. Practum® has everything you need for standard applications.

- Cost effective performance
- Supervisor lock prevents changes being made to settings
- Plug and play technology eases the transfer of data from your balance to a printer or a computer
- Includes external calibration weight
- Pan dims: Ø 90mm



Code	Alt Ref	Capacity, g	Readability, mg	Price
<b>BAL7084</b>	PRACTUM124-1S	120	0.1	£1710.00
<b>BAL7086</b>	PRACTUM224-1S	220	0.1	£1810.00