



## Cytiva Blotting Papers



High purity, smooth blotting paper which exhibits high and uniform absorbency. Made of pure cellulose produced entirely from the highest quality cotton linters with no additives of any kind to ensure that no contamination will occur during the transfer steps.

- Grade GB003: a general purpose blotting paper (0.8mm) made from pure raw materials with a high absorbency used as a membrane gel support, recommended for the lysis/denaturation of colony or plaque lifts and Western blots
- Grade GB005: a thick (1.5mm) highly absorbent paper recommended for applications where fewer layers of blotting paper must still ensure a high capacity, recommended for semi-dry blotting of proteins

Code	Alt Ref	Grade	Dims, w x l, mm	Pack	Price
<b>MOL3107</b>	10426880	GB003	100 x 100	50	<b>£24.00</b>
<b>MOL3117</b>	10426890	GB003	300 x 600	25	<b>£86.00</b>
<b>MOL3118</b>	10426972	GB005	150 x 150	25	<b>£69.00</b>
<b>MOL3119</b>	10426994	GB005	580 x 580	25	<b>£346.00</b>



Quality, value  
laboratory  
equipment for  
everyday use



## Research Standard Products



Tailored range of high  
quality equipment,  
consumables and  
reagents for research