



Sartorius Entris® II Essential Line Analytical Balances

SARTORIUS

Every Sartorius balance offers quality, value and consistency. No matter what you're weighing, the new Entris® II is always the right choice. Offering unrivalled value at a budget price, this is the only balance in its class featuring isoCAL, LED touch technology and 12 built in applications. Backed by almost 150 years of German engineering expertise, and available in 40 different models, you'll easily find an Entris® II balance that exactly meets your specific weighing needs.

- Highly accurate results year in year out – guaranteed via Monolithic weigh cell technology invented by Sartorius
- Fastest stabilisation time in its class using state of the art weighing sensors
- Optional isoCAL (internal calibration and adjustment) for total assurance of accurate weighing results
- Hybrid screen for excellent readability and use – combines intuitive, wear resistant LED and touch technology
- 12 built in applications with GLP/GMP compliant printout/data output
- USB Type C and RS232 9 pin interfaces
- Password protection for secure operation
- Pan dims: Ø 90mm

Code	Alt Ref	Capacity, g	Readability, mg	Calibration	Price
BAL1274	BCE64-1S	60	0.1	External	£1470.00
BAL1208	BCE64i-1S	60	0.1	Internal	£1712.00
BAL1276	BCE124-1S	120	0.1	External	£1560.00
BAL1218	BCE124i-1S	120	0.1	Internal	£1802.00
BAL1278	BCE224-1S	220	0.1	External	£1670.00
BAL1272	BCE224i-1S	220	0.1	Internal	£1892.00



Sartorius Entris® II Advanced Line Analytical Balances

SARTORIUS

No matter what you are weighing, the Entris® II Advanced balance range is always the right choice, offering unrivalled value and backed by almost 150 years of German engineering expertise. Featuring monolithic weigh cell technology for highly accurate results year in and year out, the fastest stabilisation time in its class using state of the art weighing sensors, a graphic touch display and 13 built in application programs.

- Built in real time level support
- Integrated protection systems with 3 configurable levels to determine valid weighing data and ensure only valid data is transferred to external devices
- High chemical resistance
- Hybrid screen for excellent readability and use combines intuitive, wear resistant LED and touch technology
- Controlled access to balance settings with user management
- Gapless documentation thanks to CalAuditTrail
- Enables multiple connection to external devices via two USB Type C interfaces
- Allows individualised printouts for up to 6 unique profiles
- Password protection for secure operation
- Verified balances are conformity assessed with EU type examination certificate
- Pan dims: Ø 90mm

Code	Alt Ref	Capacity, g	Readability, mg	Calibration	Verified	Price
BAL4622	BCA64-1S	60	0.1	External	No	£3192.00
BAL4624	BCA64i-1S	60	0.1	Internal	No	£3400.00
BAL4678	BCA64i-1CEU	60	0.1	Internal	Yes	£3868.00
BAL4564	BCA124-1S	120	0.1	External	No	£3245.00
BAL4566	BCA124i-1S	120	0.1	Internal	No	£3500.00
BAL4642	BCA124i-1CEU	120	0.1	Internal	Yes	£3912.00
BAL4580	BCA224-1S	220	0.1	External	No	£3557.00
BAL4582	BCA224i-1S	220	0.1	Internal	No	£3868.00
BAL4660	BCA224i-1CEU	220	0.1	Internal	Yes	£4286.00
BAL4592	BCA324i-1S	320	0.1	Internal	No	£4434.00
BAL4666	BCA324i-1CEU	320	0.1	Internal	Yes	£4917.00