

## OHAUS Adventurer® Analytical Balances



Striking the ideal balance between inventive features and functional, uncomplicated weighing capabilities, the OHAUS Adventurer® incorporates all of the applications necessary for routine weighing and measurement activities. With a colour touchscreen, GLP/GMP compliance capabilities, two USB ports, and much more, Adventurer® is the most complete balance in its class.

- Fast stabilisation and reliable operation
- AutoCal™ option to ensure solid weighing performance and accurate, repeatable results for standard laboratory applications
- Wide viewing angle colour touchscreen provides access to all of Adventurer®'s applications
- Dual USB ports with GLP/GMP tracking capabilities assist in the monitoring and reporting of data
- Two piece, top mounted draftshield doors reduce the footprint of the doors in the rear of the balance when open, and open wide on both sides in order to provide ample room to place samples
- Dims, w x d x h: 354 x 230 x 340mm
- Pan dims: Ø 90mm



Code	Alt Ref	Model	Capacity, g	Readability, mg	Calibration	Price
<b>BAL0011</b>	30122657	AX124/E	120	0.1	External	£2020.00
<b>BAL0010</b>	30122610	AX124	120	0.1	Internal	£2285.00
<b>BAL0013</b>	30122612	AX224/E	220	0.1	External	£2135.00
<b>BAL0012</b>	30122611	AX224	220	0.1	Internal	£2405.00
<b>BAL0014</b>	30122613	AX324	320	0.1	Internal	£2425.00

## OHAUS Explorer® Analytical Balances



Featuring AutoCal™ automatic internal calibration whilst enhanced vibration filtering ensures fast and stable weighing results. A colour touch screen adjustable angle display with intuitive icon driven application software provides access to 14 application modes make the Explorer® easy to use whilst the separation of base and display maximise flexibility.

- High resolution multilingual colour VGA touch display with icons for simple operation
- Programmable multifunctional touch free sensors for cleanliness, speed and safety
- Over 14 application modes including pipette adjustment
- Library for storage of personal application settings
- RS232 and USB as standard, plus 3rd connection slot (ethernet or 2nd RS232)
- Full GLP/GMP including date and time
- Supplied with data transfer software as standard
- Modular display design with adjustable angle (attach to base, wall or optional tower)
- AutoCal™ automatic calibration
- Level assist feature
- Integral weigh below hook
- Automatic draftshield models available
- Pan dims: Ø 90mm



Code	Alt Ref	Model	Capacity, g	Readability, mg	Auto Door	Price
<b>BAL5100</b>	83021331	EX124	120	0.1	No	£3230.00
<b>BAL5102</b>	83021332	EX224	220	0.1	No	£3500.00
<b>BAL5106</b>	83021335	EX324	320	0.1	No	£4050.00
<b>BAL5101</b>	30061976	EX124/AD	120	0.1	Yes	£4070.00
<b>BAL5103</b>	30061977	EX224/AD	220	0.1	Yes	£4340.00
<b>BAL5105</b>	30061978	EX324/AD	320	0.1	Yes	£4890.00