

## Buchi F-100 and F-105 Recirculating Chillers Continued



Specifications	
Max Power Consumption, W	850
Heating Emission, W	700
Refrigerant	R 134a
Tank Volume, L	3
Max Pump Pressure, bar	0.6
Flow Rate, L/min	2.5
Dims, w x d x h, mm	280 x 420 x 500

Code	Alt Ref	Model	Display	Temp Range, °C	Cooling Capacity, W	Price
EVA4188	11056460	F-100	No	10 (fixed)	@10°C: 300	£1893.00
EVA4190	11060002	F-105	Digital	-10 to +25	@15°C: 530, @10°C: 390, @0°C: 120	£2487.00

## Accessories

Code	Alt Ref	Description	Price
EVA4198	037742	Distribution adapter for operating two extraction units with a single chiller	£422.00

## JULABO F/FL Series Recirculating Coolers



Requiring very little space, the efficiency of these coolers makes them an environmentally friendly and economical. The F series is a great way to replace costly tap water and ideal for basic cooling tasks, and the FL series offer increased flexibility and power for more demanding applications.

- Environmentally friendly operation with low energy consumption
- Compact design
- Splashproof membrane keypad with LED temperature display
- Straightforward filling and draining
- Filling level indicator
- May be used with water, water/glycol, JULABO Thermal G
- Temperature stability:  $\pm 5^{\circ}\text{C}$
- 1PLUS warranty offers an extension to the normal warranty of up to 24 months, limited to a maximum of 10000 working hours



Code	Alt Ref	Model	Temp Range, °C	Cooling Capacity, kW	Flow Rate, L/min	Pump Pressure, bar	Price
BAT8050	9620025	F250	-10 to 40	0.25	15	0.35	£2580.00
BAT8055	9620050	F500	0 to 40	0.5	24	0.5	£3674.00
BAT8060	9620100	F1000	0 to 40	1	23	1	£4595.00
BAT2500	9660003	FL300	-20 to 40	0.3	15	0.35	£3659.00
BAT2502	9661006	FL601	-20 to 40	0.6	23	1	£4921.00
BAT2504	9661012	FL1201	-20 to 40	1.2	23	1	£6249.00