

Cole-Parmer® EME/EMEA Multi-Position Electromantles Continued

Cole-Parmer®

Code	Alt Ref	Model	No. of Places	Flask Size, mL	Stirrer	Power, W	Dims, w x d x h, mm	Price
EXT1000	03076-01	EME30100/CEB	3	100	No	180	640 x 270 x 110	£864.00
EXT1001	03076-13	EME30100/CE	3	100	Yes	220	640 x 270 x 110	£1585.00
EXT1002	03076-04	EME30250/CEB	3	250	No	450	640 x 270 x 110	£973.00
EXT1003	03076-16	EMEA30250/CE	3	250	Yes	490	640 x 270 x 110	£1585.00
EXT1004	03076-07	EME30500/CEB	3	500	No	600	755 x 285 x 130	£1114.00
EXT1005	03076-19	EMEA30500/CE	3	500	Yes	640	755 x 285 x 130	£1702.00
EXT1006	03076-10	EME31000/CEB	3	1000	No	900	755 x 285 x 130	£1108.00
EXT1007	03076-22	EMEA31000/CE	3	1000	Yes	940	755 x 285 x 130	£1742.00
EXT1010	03076-25	EME60100/CEB	6	100	No	420	1210 x 270 x 110	£1375.00
EXT1011	03076-37	EMEA60100/CE	6	100	Yes	500	1210 x 270 x 110	£2688.00
EXT1012	03076-28	EME60250/CEB	6	250	No	900	1210 x 270 x 110	£1583.00
EXT1013	03076-40	EMEA60250/CE	6	250	Yes	980	1210 x 270 x 110	£2484.00
EXT1014	03076-31	EME60500/CEB	6	500	No	1200	1450 x 285 x 130	£1725.00
EXT1015	03076-43	EMEA60500/CE	6	500	Yes	1280	1450 x 285 x 130	£2999.00
EXT1017	03076-46	EMEA61000/CE	6	1000	Yes	1880	1450 x 285 x 130	£3086.00

Quickfit® Soxhlet Extractors



- Conforms to BS 2071
- For extracting soluble components or impurities from solids via repeated soaking/siphoning with solvent
- Accurately formed concentric siphon tube for reliable performance



Code	Cone Size	Socket Size	Recommended Thimble Size, mm	Nominal Capacity (Siphoning Volume), mL	Price
QEX555100	34/35	34/35	28 x 120	100	£91.00
QEX563	24/29	40/38	33 x 100	100	£88.00
QEX583	24/29	50/42	41 x 123	200	£109.00

Visit our dedicated safety subsite

Bringing you easy to navigate categorisation, extensive choice and live pricing. Visit our safety subsite to view our ever increasing range of safety products.

scientificlabs.co.uk/safety (UK & Northern Ireland)

scientificlabs.ie/safety (Republic of Ireland)

scientificlabs.com/safety (International)

