

Heathrow Scientific Parafilm® M Sealing Film Dispensers



- Store, dispense, and cut Parafilm® M, tape, labels and Tough-Tags®
- Holds 2 x 50mm rolls, or 1 x 100mm roll
- Built in safety blade slides and cuts with no danger to fingertips, eliminates the need for scissors
- Translucent dispenser allows users to track when film is getting low, and replace
- Storage bin on back of the dispenser keeps pens and other small accessories handy
- Translucent colours help prevent film from yellowing
- Dims, w x d x h: 120 x 156 x 171mm

Code	Alt Ref	Description	Pack	Price
FIL1034	HS234525B	Blue	Each	£76.00
FIL1038	HS234525C	Green	Each	£58.00
FIL1028	HS234525D	Red	Each	£58.00
FIL1035	HS234525R	Replacement blades	5	£17.20



Bringing you the established brands across the protein research workflow


Sera-Mag™ nucleic acid sample and library preparation

HyClone™ cell culture media and sera

Whatman™ laboratory filtration products

ÄKTA™ chromatography columns, resins, and systems

Amersham™ systems, membranes, films, and reagents



cytiva.com

Cytiva and the Drop logo are trademarks of Life Sciences IP Holdings Corp. or an affiliate doing business as Cytiva, ÄKTA, Amersham, HyClone, Sera-Mag, and Whatman are trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.

© 2022 Cytiva

For local office contact information, visit cytiva.com/contact

CY21794-07Mar22-AD



Read
more
here