

SLS Select Qualitative Filter Paper Discs



Grade A

The most widely used general purpose filter paper for general applications in laboratory analysis and liquid clarification. Medium filtration velocity and retention. Whatman Grade 1 equivalent. 160µm thickness, 5 - 13µm nominal retention.

Supplied in packs of 100

Code	Ø, mm	Price
SLS2028	55	£4.20
SLS2029	70	£4.60
SLS2031	90	£5.30
SLS2032	110	£6.20
SLS2033	125	£6.40
SLS2034	150	£9.10
SLS2035	185	£13.60
SLS2036	240	£21.90



Grade B

Qualitative filter paper, suitable for the fast filtration of liquids with fine particles. Whatman Grade 2 equivalent. 160µm thickness, 5 - 8µm nominal retention.

Supplied in packs of 100

Code	Ø, mm	Price
SLS2054	90	£7.50
SLS2055	110	£8.90
SLS2056	125	£11.20
SLS2057	150	£14.80



Grade D

Qualitative filter paper, fast filtering, suitable for coarse precipitates such as hydroxides. Very useful as a rapid filter for routine clean-up of biological fluids or organic extracts. Whatman Grade 4 equivalent. 200µm thickness, 12 - 15µm nominal retention.

Supplied in packs of 100

Code	Ø, mm	Price
SLS2088	90	£7.00
SLS2089	110	£9.00
SLS2090	125	£11.10
SLS2091	150	£15.10



Grade F

Qualitative filter paper, suitable for the filtration of fine precipitates such as nickel sulphide, lead dioxide and calcium fluoride. Whatman Grade 6 equivalent. 140µm thickness, 3 - 5µm nominal retention.

Supplied in packs of 100

Code	Ø, mm	Price
SLS2112	90	£7.90
SLS2113	110	£9.90
SLS2114	125	£12.50
SLS2115	150	£16.60

