

Sartorius Vivaspin® Centrifugal Concentrators

SARTORIUS

With their broad range of volume capacities and MWCOs, Vivaspin® centrifugal concentrators offer the optimal solution to any concentration or buffer exchange application. Due to their particular vertical membrane design they achieve highest flow rates and recoveries. The patented dead-stop technology allowing direct recovery of the concentrate without further re-spin assures incomparably safe sample handling and convenience.

- Vertical, large membrane design
- Integrated dead-stop-volume
- Low-binding materials
- No re-spin necessary
- Broad volume range and wide MWCO range
- Fast, minimal membrane blocking
- No risk (to concentrate to dryness)
- High recovery
- Easy handling, fast



Vivaspin® 500, PES Membrane

Volume: 0.1 - 0.5mL.

Code	Alt Ref	MWCO, KDa	Pack	Price
FIL8546	VS0191	3	25	£111.00
FIL8548	VS0192	3	100	£370.00
FIL8550	VS0111	5	25	£111.00
FIL8552	VS0112	5	100	£370.00
FIL8554	VS0101	10	25	£126.00
FIL8556	VS0102	10	100	£370.00
FIL8558	VS0121	30	25	£111.00
FIL8560	VS0122	30	100	£370.00
FIL8562	VS0131	50	25	£111.00
FIL8564	VS0132	50	100	£370.00
FIL8566	VS0141	100	25	£111.00
FIL8568	VS0142	100	100	£370.00
FIL8569	VS0151	300	25	£112.00
FIL8571	VS0152	300	100	£370.00
FIL8573	VS0161	1000	25	£111.00
FIL8575	VS0162	1000	100	£370.00
FIL8577	VS0171	0.2µm	25	£111.00
FIL8579	VS0172	0.2µm	100	£368.00



Vivaspin® 2, PES Membrane

Volume: 0.4 - 2mL.

Code	Alt Ref	MWCO, KDa	Pack	Price
FIL8597	VS0291	3	25	£145.00
FIL8599	VS0292	3	100	£487.00
FIL8600	VS0211	5	25	£145.00
FIL8602	VS0212	5	100	£487.00
FIL8604	VS0201	10	25	£183.00
FIL8606	VS0202	10	100	£409.00
FIL8608	VS0221	30	25	£145.00
FIL8610	VS0222	30	100	£487.00
FIL8612	VS0231	50	25	£145.00
FIL8614	VS0232	50	100	£487.00
FIL8616	VS0241	100	25	£145.00
FIL8618	VS0242	100	100	£487.00
FIL8652	VS0251	300	25	£145.00
FIL8621	VS0252	300	100	£487.00
FIL8623	VS0261	1000	25	£145.00
FIL8625	VS0262	1000	100	£487.00
FIL8627	VS0271	0.2µm	25	£145.00
FIL8629	VS0272	0.2µm	100	£555.00