



## Whatman™ 3 Piece Filter Funnels

cytiva  
Whatman

- Ideal for quick and easy filtration using high-efficiency glass microfibre filters
- Consists of a glass funnel and base, acrylic support platform, and 3 easy-to-apply metal clips.
- Range of diameters, funnel sizes, reservoir volumes and alternate support platforms provides a flexible system to meet the varied and unique needs of busy laboratories
- Chemically resistant construction of glass cylinder and base with a polymeric membrane support
- Provided with acrylic support and optional polypropylene or PTFE supports available separately for some sizes
- Clamps at edges to prevent peripheral loss due to fluid bypass around the edge of the filter
- Simplicity of design facilitates easy and thorough cleaning between use

Code	Alt Ref	Reservoir Capacity, mL	Filter Paper Size, mm	Price
<b>FUN3100</b>	1950-002	16	25	£165.00
<b>FUN3102</b>	1950-004	36	47	£189.00
<b>FUN3104</b>	1950-007	115	70	£228.00
<b>FUN3106</b>	1950-017	210	70	£251.00
<b>FUN3108</b>	1950-027	400	70	£393.00
<b>FUN3110</b>	1950-009	200	90	£262.00
<b>FUN3112</b>	1950-012	530	125	£533.00

## Alternative Support Plates

PTFE.

Code	Alt Ref	Filter Paper Size, mm	Price
<b>FUN3126</b>	1950114	47	£71.00
<b>FUN3130</b>	1950117	70	£90.00



## SLS Select Syringe Filters



The highly asymmetrical pore structure of PES membranes offers an excellent loading capacity and high flow rate for our syringe filters which are naturally hydrophilic.

- Chemical stability: pH 4 – 8
- Recommended for: biology, preparation of protein samples, ultra-fast filtration of tissue culture media
- Diameter: 25mm
- Sterile

**Supplied in packs of 50 (individually wrapped)**

Code	Pore Size, µm	Price
<b>SLS3011</b>	0.22	£67.00
<b>SLS3012</b>	0.45	£67.00