

Sartorius Minisart® RC Syringe Filters

SARTORIUS

Non sterile Minisart® regenerated cellulose (RC) syringe filters are supplied in a pouch in one pack. This reduces waste and makes the handling easier whenever the application allows clarification under non sterile conditions, like this is the case for the majority of sample preparation steps prior to HPLC or GC analysis. The RC membrane is the best option for reliable and leachable free filtration. The very low non specific binding characteristics of this membrane makes it a perfect choice for aqueous samples and solvents. In addition, unlike PVDF filters, those with RC membranes are suitable for DMSO and other amides, or for ketones, esters, and ethers. The housing of the filter is made of polypropylene (PP) which is compatibility with harsh solvents, acids and bases.

- A pouch in one pack for less waste
- Minimum extractables or leachables
- Very low non specific binding
- Outperforming chemical compatibility



Code	Alt Ref	Membrane Ø, mm	Pore Size, µm	Connector Outlet	Colour	Sterility	Pack	Price
FIL6652	17821K	4	0.2	Male luer slip	Blue tray	Non sterile	50	£75.00
FIL6654	17821Q	4	0.2	Male luer slip	Blue tray	Non sterile	500	£633.00
FIL6656	17822K	4	0.45	Male luer slip	Yellow tray	Non sterile	50	£74.00
FIL6658	17822Q	4	0.45	Male luer slip	Yellow tray	Non sterile	500	£633.00
FIL6660	17761K	15	0.2	Male luer slip	White, printed	Non sterile	50	£105.00
FIL6662	17761Q	15	0.2	Male luer slip	White, printed	Non sterile	500	£723.00
FIL6668	17764K	25	0.2	Male luer slip	White, printed	Non sterile	50	£117.00
FIL6670	17764Q	25	0.2	Male luer slip	White, printed	Non sterile	500	£781.00
FIL6672	17765K	25	0.45	Male luer slip	White, printed	Non sterile	50	£117.00
FIL6674	17765Q	25	0.45	Male luer slip	White, printed	Non sterile	500	£781.00

Sartorius Minisart® HY Syringe Filters

SARTORIUS

Bacterial cell removal or particulate removal from liquids, including medical drugs, can easily be performed. Minisart® HY with hydrophobic PTFE have the CE mark of conformity and are frequently used for sterile filtration of aqueous and oily ear or eye solutions and other drugs.

- 6mm, PTFE membrane
- CE marked air and gas filtration and venting
- Membrane diameter: 26mm



Code	Alt Ref	Pore Size, µm	Connector Outlet	Colour	Sterility	Pack	Price
FIL6720	16596HYK	0.2	Male luer lock	Clear	Ethylene oxide	50	£131.00
FIL6722	16596HYQ	0.2	Male luer lock	Clear	Non sterile	500	£851.00