

Whatman™ GD/X® Syringe Filters



Specifically designed for high particulate loaded samples. Contains a prefiltration stack to eliminate sample contamination. Enables high loading capacity, fast flow rates and eliminates back pressure allowing difficult samples to be filtered with less hand pressure. GD/X® syringe filters can process three to seven times more sample volume than unprotected membranes. Especially useful for particulate-laden samples found in: dissolution testing, content uniformity, concentration analysis, food and environmental samples and composite assay.

- 13mm filters for samples up to 10mL
- 25mm filters for samples ≥10mL
- Sterile options
- Pigment-free polypropylene housing
- Novel prefiltration stack of GMF 150 (graded density) and GF/F glass
- Microfibre media



13mm Syringe Filters

Code	Alt Ref	Membrane	Sterility	Pore Size, µm	Pack	Price
FIL2541	6870-1302	Nylon	Non sterile	0.2	150	£389.00
FIL2543	6870-1304	Nylon	Non sterile	0.45	150	£389.00
FIL2548	6871-1304	Nylon	Non sterile	0.45	1500	£2354.00
FIL2547	6872-1304	PVDF	Non sterile	0.45	150	£389.00
FIL2550	6873-1304	PVDF	Non sterile	0.45	1500	£2386.00
FIL2549	6874-1302	PTFE	Non sterile	0.2	150	£395.00
FIL2552	6875-1302	PTFE	Non sterile	0.2	1500	£2386.00
FIL2551	6874-1304	PTFE	Non sterile	0.45	150	£389.00
FIL2554	6875-1304	PTFE	Non sterile	0.45	1500	£2354.00
FIL2553	6876-1302	PES	Non sterile	0.2	150	£393.00
FIL2555	6878-1302	PP	Non sterile	0.2	150	£393.00
FIL2557	6878-1304	PP	Non sterile	0.45	150	£393.00
FIL2559	6880-1304	CA	Non sterile	0.45	150	£389.00
FIL2560	6882-1316	GMF GF/A*	Non sterile	1.6	150	£389.00
FIL2562	6884-1310	GMF GF/B*	Non sterile	1.0	150	£389.00
FIL2564	6886-1312	GMF GF/C*	Non sterile	1.2	150	£392.00
FIL2566	6888-1327	GMF GF/D*	Non sterile	2.7	150	£389.00

* Particle retention rating

25mm Syringe Filters

Code	Alt Ref	Membrane	Sterility	Pore Size, µm	Pack	Price
FIL2570	68692502	Nylon	Non sterile	0.2	150	£416.00
FIL4150	68702502	Nylon	Non sterile	0.2	150	£411.00
FIL2572	68692504	Nylon	Non sterile	0.45	150	£413.00
FIL4152	68702504	Nylon	Non sterile	0.45	150	£411.00
FIL2574	68712550	Nylon	Non sterile	5.0	1500	£2442.00
FIL4154	68702550	Nylon	Non sterile	5.0	150	£411.00
FIL2590	68762502	PES	Non sterile	0.2	150	£411.00
FIL2598	68972502	PES	Sterile	0.2	500	£1925.00
FIL2594	69052502	PES	Non sterile	0.2	1500	£2442.00
FIL2596	69052504	PES	Non sterile	0.45	1500	£2428.00
FIL4166	68762504	PES	Non sterile	0.45	150	£415.00
FIL4156	68722502	PVDF	Non sterile	0.2	150	£411.00
FIL4158	68722504	PVDF	Non sterile	0.45	150	£411.00
FIL2576	68732502	PVDF	Non sterile	0.45	1500	£2428.00
FIL4160	68742502	PTFE	Non sterile	0.2	150	£411.00